


GPC / SEC / HPLC Columns

Made by AppliChrom



Product Catalog 2021

Sugars

Degradation Products

Polymers

Biopolymers

Peak-Recycling-GPC



Table of Contents

I.	About Us	1
II.	AppliChrom ABOA DMAC-Phil-P	3
III.	AppliChrom ABOA DMSO-Phil-P	5
IV.	AppliChrom ABOA SuperOH-P	12
V.	AppliChrom ABOA SugarSep	20
VI.	AppliChrom ABOA SugarSep-Oligo	24
VII.	AppliChrom ABOA StyDiViBe-P	26
VIII.	AppliChrom ABOA CatPhil-P	34
IX.	AppliChrom ABOA CatPhil-P-JLJ	37
X.	AppliChrom HFIP-Phil-P	39
XI..	AppliChrom Aceton-Aq-P	42
XII.	AppliChrom Peak-Recycling-GPC	44
XIII.	AppliChrom ABOA Series GPC/SEC Column Guide	46
XIV.	AppliChrom HPLC / GPC Column User Guide	48
XV.	AppliChrom General Terms and Conditions	49
XVI.	AppliChrom Application List of Application / Product Catalog 2021	50



About Us

AppliChrom team and customers grow together by supporting chromatographers in actual challenging and often difficult separations. You can count on AppliChrom with its own production site for chromatography media, columns and applications and more than 50 years of HPLC/GPC experience at the traditional historical place where 1850 chromatography was first invented by Professor Runge (Dipl.-Ing. (FH) Susanne & Dr. Christian Dauwe).

13 Reasons for AppliChrom

1) We take it personally

Especially the support of our customers. From the analytical question up to the chromatographic process we discuss and accompany you on the way to your chromatographic result.

2) Application and Technical Support from one Source

As your chromatography partner we organize and look at the whole chromatographic process – also under optimal use of the HPLC/GPC/SEC-devices already in use.

3) Increasing Quality

You are looking for the best way to get your optimum chromatographic result in HPLC, GPC or SEC. We do not only offer high quality and highly reproducible products for high resolution and lifetime, but also a suitable consultation to your individual chromatographic task.

4) Saving Time

Chromatographers who start developing a method have to finally produce a reliable method that give complete information about all components of interest from a sample. Maybe the method has to be compatible to different kinds of chromatography systems or environments later – taking these aspects into account from the beginning saves much time and cost for our customers. Our product specialists are glad to assist you from the beginning of a project to make you successful.

5) Experience

AppliChrom's experience of own production of chromatography media, columns and applications (HPLC, GPC and SEC) allows us to support you with the background knowledge of the chromatography media producer with staff having more than 50 years of relevant experience.

6) Individual

From the beginning of your request we integrate your individual wishes, ideas and needs into the support to your optimal chromatographic solution. Resolution, selectivity, compatibility to existing systems, time saving or all of this – please ask us. AppliChrom has own application laboratories, own manufacturing capacities and own development unions to support you with your individual needs.

7) Innovative

AppliChrom chromatography media are typically based on special customers' needs that were covered via the AppliChrom development process unit finishing with the process validation to high end products that are available now for your reliable product solutions. Thus many very unique and powerful solutions are offered now for HPLC, SEC and GPC by AppliChrom.

8) Reproducibility

AppliChrom works at the headquarter Oranienburg has its own quality control laboratories, the production of chromatography columns in small and large series but also units to produce chromatography media in small scale and for larger bulk amounts. The quality of each batch and each column is controlled to ensure the highest level of reproducibility, column for column, batch to batch and year after year.

9) OEM business

AppliChrom supplies well-chosen OEM customers with individual chromatography bulk or column lines. The production also contains customized lines for special needs. You do not find these products necessarily in our public catalogue. If you have a special need – please speak with us.

10) Chromatography – Origin from Oranienburg, innovations from Oranienburg

Maybe it is interesting for you when and where chromatography was invented. It was in 1850 in Oranienburg where professor Runge invented, published and commercialized chromatography products first in the world – it was an early form of paper chromatography – he called it: „Bilder die sich selber malen“ known also as Runge pictures. These are the roots of chromatography and the heritage of Oranienburg that we develop here at Oranienburg to innovative high end products to support you with your specific chromatographic solution.

11) Made in Europe

AppliChrom brand products are quality products made in Germany, coming from Oranienburg located in the Berlin region. Surrounded by several universities and institutes dealing with life and material science we get a continuous input of actual scientific challenges to get inspirations for needed new high end products. Anyway we are glad to hear what AppliChrom can do for you.

12) Carefree

Chromatographic separation should not be an adventure for the customer. We serve you from your request with the needed individual assistance in the beginning up to the after sales service/consultation.

13) Satisfaction

At the end AppliChrom's goal is to only have satisfied customers. This is our endeavour.

AppliChrom's core competence is the development and production of high quality liquid chromatography media in analytical and semipreparative scale. For this reason AppliChrom is using different reactor lines for producing several unique kinds of chromatography media. Especially for chromatography (HPLC, SEC, GPC) of biopolymers, oligomers, degradation products, fermentation products AppliChrom serves its customers with a unique product line.

AppliChrom branded products are quality products „Made in Germany“, these products are developed and produced in Germany and are submitted to a strict high class quality control.

The headquarter of AppliChrom is located in Oranienburg near Berlin. Oranienburg is the town where Friedlieb Ferdinand Runge invented the early form of paper chromatography in 1850.

Analytical chromatography is a process of separation, quantification and identification of chemical or biological compounds. Chromatography is used e.g. for research and analysis of food, drinks, pharmaceuticals, metabolomic processes, material sciences, environment and chemicals.

With our more than 5.000 products we supply chromatography columns (HPLC, SEC, GPC) for all current applications. Therefore, numerous institutes, universities, manufacturers, as well as our other customers appreciate the high quality, lifetime, reproducibility and reliability of AppliChrom products.

We extend our program steadily with new own chromatographic developments which fulfill the special demands of our customers to secure an economic workflow at the customers' site.

We are also glad about your inquiry!

Important facts about AppliChrom

Founded: 2006, first market presence in 2009

Company founder: Dipl.-Ing. (FH) Susanne Dauwe (Chemical engineer)

Research, development and chromatography applications: Dr. Christian Dauwe (Chemist)

Products:


Silica- and polymeric based liquid chromatography media (HPLC, GPC, SEC)

HPLC-columns – C8, C18, HILIC & hybrid, ion exchanger and ion exclusion columns and media

Special easy to use HPLC columns for sugars, alcohols, acids (HILIC, ion exclusion, ion exchangers)

SEC/GPC columns for water soluble biopolymers

SEC/GPC columns for organic soluble polymers



Monument of Dr. Runge in front of the former Laboratory in Oranienburg

AppliChrom ABOA DMAc-Phil-P

AppliChrom GPC columns for GPC analyzes organic (DMAc, DMF and NMP)

Spherical porous polymeric GPC chromatography phases


- wide range of pore and particle sizes
- optimal for 0.5ml / min flow rate with 8mm ID columns (40-50°C)
- Pressure stability 50-150bar, depending on the porosity
- high resolution due to high pore volume
- long service life
- high reproducibility
- high purity of the GPC column for good interaction-free GPC
- Please avoid: drying out
- Molar mass range: 100-1,000,000Da
- Examples: PMMA, PAN, cellulose, DMF and DMAc soluble polymers

AppliChrom ABOA DMAc-Phil-P molar mass range:


P-100	100-2.500Da
P-200	100-20.000Da
P-250	100-70.000Da
P-250JLD*)	100-100.000Da
P-300	1000-300.000Da
P-350	1.500-1.000.000Da

*) JLD Pore: extended linear range

To cover a very wide range of molecular sizes, GPC columns of appropriate porosity can be combined.



GPC Calibration Curves AppliChrom ABOA DMAc-Phil



Long linear calibration range

Analyte:

PMMA Standards

Column:

AppliChrom ABOA DMAcPhil-P-**100**
 AppliChrom ABOA DMAcPhil-P-**200**
 AppliChrom ABOA DMAcPhil-P-**250**
 AppliChrom ABOA DMAcPhil-P-**250-JLD**
 AppliChrom ABOA DMAcPhil-P-**300**
 AppliChrom ABOA DMAcPhil-P-**350**

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

H₂O

Flow:

0,5ml/min

Temperature:

45°C

Detection:


RI

Injection:

20µl sample

AppliChrom ABOA DMAC-Phil-P

Catalog #	Description	Dimension	Separation Range
SADCP1002508 SADCP1003008 SADCP100508 SADCP100308	AppliChrom ABOA DMAC-Phil-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SADCP2002508 SADCP2003008 SADCP200508 SADCP200308	AppliChrom ABOA DMAC-Phil-P-200	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-20.000Da
SADCP2502508 SADCP2503008 SADCP250508 SADCP250308	AppliChrom ABOA DMAC-Phil-P-250	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-70.000Da
SADCP2502508J SADCP2503008J SADCP250508J SADCP250308J	AppliChrom ABOA DMAC-Phil-P-250	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-100.000Da
SADCP3002508 SADCP3003008 SADCP300508 SADCP300308	AppliChrom ABOA DMAC-Phil-P-300	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-300kDa
SADCP3502508 SADCP3503008 SADCP350508 SADCP350308	AppliChrom ABOA DMAC-Phil-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-1.MioDa
SADCP4002508 SADCP4003008 SADCP400508 SADCP400308	AppliChrom ABOA DMAC-Phil-P-400	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10kDa-5 MioDa
SADCP4502508 SADCP4503008 SADCP450508 SADCP450308	AppliChrom ABOA DMAC-Phil-P-450	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100 kDa→10MioDa



Reproducibility

AppliChrom ABOA DMSO-Phil-P

Special GPC/SEC media for fast, accurate and robust GPC-analysis in DMSO


For GPC / SEC analysis in DMSO, examples:

- amylose, amylopectin, starch
- urea-formaldehyd resins (UF-resins)
- melamin-urea-formaldehyd resins (MUF-resins)
- lignins, humic substances, humic acids, coniferous wood bark essences
- polysaccharide, polysaccharid derivatives
- poly(N-isopropylacrylamid) PNIPA
- poly-vinylpyridin
- calibration: pullulan, dextran, polyvinylpyridin et al.


Advantages

- optimized for DMSO-GPC applications
- interactionfree pure GPC/SEC
- easy, reliable and robust GPC/SEC-calibration by dextrans, pullulans et al.
- low column bleeding for low detectornoise for improved lightscattering or viscosity detection
- 12 μ particle technology for low backpressure
- large porevolume and optimized mass transfer for polymers giving optimized resolution
- low costs caused by long livetime of column – result of combination of optimized proprietary particle and packingtechnology.

Measuring range of selected AppliChrom ABOA DMSO-Phil-P GPC/SEC-columns, 300x8mm, poresizes 100, 200, 250, 350, 400




Different pore sizes available



Analyte: Lignin conifer bark extracts

Column: AppliChrom ABOA DMSO-Phil-P-250
AppliChrom ABOA DMSO-Phil-P-350


Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,5ml/min
Temperature: 80°C
Detection: RI



Analyte: Spruce bark extract

Column: AppliChrom ABOA DMSO-Phil-P-250

Dimension: 2 x 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,5ml/min
Temperature: 80°C
Detection: RI




Analyte: Polysaccharide

(M ca. 70kDa)

Column: AppliChrom ABOA DMSO-Phil-P-250
AppliChrom ABOA DMSO-Phil-P-350

Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,5ml/min
Temperature: 80°C
Detection: RI




Analyte: Dextran 650

Dextran from Leuconostoc spp.,
M = 450.000-650.000Da + fructose


Column: AppliChrom ABOA DMSO-Phil-P-200
AppliChrom ABOA DMSO-Phil-P-250
AppliChrom ABOA DMSO-Phil-P-350

Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,5ml/min
Temperature: 80°C
Detection: RI

**Analyte: Pea starch**

Column: AppliChrom ABOA DMSO-Phil-P-200
AppliChrom ABOA DMSO-Phil-P-250
AppliChrom ABOA DMSO-Phil-P-350


Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,5ml/min
Temperature: 80°C
Detection: RI

**Analyte: MUF-resin**

completely DMSO-soluble melamin-urea formaldehyd resin (MUF-resin)

Column: AppliChrom ABOA DMSO-Phil-P-200
AppliChrom ABOA DMSO-Phil-P-250
AppliChrom ABOA DMSO-Phil-P-350


Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,5ml/min
Temperature: 80°C
Detection: RI
Injection: 50µl sample

**Analyte: UF-resin**

DMSO-soluble urea formaldehyd resin

Column: AppliChrom ABOA DMSO-Phil-P-200
AppliChrom ABOA DMSO-Phil-P-350


Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,5ml/min
Temperature: 80°C
Detection: RI

**Analyte: Poly(N-isopropylacrylamide)**

Further denominations: PNIPA,
PNIPAAm, NIPA, PNIPAA or PNIPAm.
CAS [25189-55-3], formula: (C₆H₁₁NO)_n
3 different PNIPA fractions

Column: AppliChrom ABOA DMSO-Phil-P-300

Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,5ml/min
Temperature: 80°C
Detection: RI

**Analyte:****Polyvinylpyridine degraded**

Further denominations: CAS 25014-15-7, $(\text{C}_7\text{H}_7\text{N})_n$
low molecular weight (oligomeric)
polyvinylpyridin fraction

Column:

AppliChrom ABOA DMSO-Phil-P-250

Dimension:

ea. 300mm x 8mm

Mobil Phase:

DMSO / 0,075M NaNO_3

Flow:


0,4ml/min

Temperature:

50°C

Detection:

RI

**Analyte:****Polyvinylpyridine fractions**

Further denominations: CAS 25014-15-7, $(\text{C}_7\text{H}_7\text{N})_n$
75.7kDa, 20.9kDa, 3.2kDa

Superposition of 3 different polyvinylpyridin fractions

Column:

AppliChrom ABOA DMSO-Phil-P-250

Dimension:

ea. 300mm x 8mm

Mobil Phase:

DMSO / 0,075M NaNO_3

Flow:


0,4ml/min

Temperature:

50°C

Detection:

RI

**Analyte:****Polyvinylalcohol M=22kDa**

Column:

AppliChrom ABOA DMSO-Phil-P-300

Dimension:

2x. 300mm x 8mm

Mobil Phase:

DMSO / 0,075M NaNO_3

Flow:


0,4ml/min

Temperature:

50°C

Detection:

RI

**Analyte:****Poly[2-(butenyl)2-oxazoline-co-ethylenimine]****M = 50.000Da**

Column:

AppliChrom ABOA DMSO-Phil-P-300

Dimension:

ea. 300mm x 8mm

Mobil Phase:

DMSO / 0,075M NaNO_3

Flow:


0,4ml/min

Temperature:

50°C


Detection:

RI




Analyte: Polybutyl methacrylate/
Poly(dimethylamino-
ethylmethacrylate) M=22kDa

Column: AppliChrom ABOA DMSO-Phil-P-300
Dimension: 2x. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,4ml/min
Temperature: 50°C
Detection: RI




Analyte: Pea protein GPC
covering the calibration range of 100-
1.000.000Da (based on dextrans)

Column: AppliChrom ABOA DMSO-Phil-P-Multipore
Dimension: 3x. 300mm x 8mm
Mobil Phase: DMSO
Flow: 0,4ml/min
Temperature: 55°C
Detection: UV 280nm
Calibration: vs. Dextran, pullulan or protein




Analyte: Soy protein GPC
covering the calibration range of
100-1.000.000Da (based on dextrans)

Column: AppliChrom ABOA DMSO-Phil-P-Multipore
Dimension: 3x. 300mm x 8mm
Mobil Phase: DMSO
Flow: 0,4ml/min
Temperature: 55°C
Detection: UV 280nm
Calibration: vs. Dextran, pullulan or protein



Analyte: Manuka honey protein GPC
covering the calibration range of
100-1.000.000Da (based on dextrans)

Column: AppliChrom ABOA DMSO-Phil-P-Multipore
Dimension: 3x. 300mm x 8mm
Mobil Phase: DMSO
Flow: 0,4ml/min
Temperature: 40°C
Detection: UV 270nm
Calibration: vs. Dextran, pullulan or protein




Analyte: Caramel color GPC

(range: 100Da-1.500.000Da)

SEC / GPC comparison of 2 Caramel
Coulor

Column: AppliChrom ABOA DMSO-Phil-P-100
AppliChrom ABOA DMSO-Phil-P-350

Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,3ml/min
Temperature: 60°C
Detection: UV 280nm




Analyte: Humic acids / humates GPC

(range: 100-1.500.000Da)

SEC / GPC comparison GPC / SEC
comparison of 2 different humic acids /
humates

Column: AppliChrom ABOA DMSO-Phil-P-Multipore

Dimension: 3x 300mm x 8mm
Mobile Phase: DMSO
Flow: 0,4ml/min
Temperature: 70°C
Detection: UV 280nm




Analyte: Melamin GPC

(range: 100-1.500.000Da)


Column: AppliChrom ABOA DMSO-Phil-P-100
AppliChrom ABOA DMSO-Phil-P-350

Dimension: ea. 300mm x 8mm
Mobil Phase: DMSO / 0,075M NaNO₃
Flow: 0,4ml/min
Temperature: 60°C
Detection: RI



AppliChrom ABOA DMSO-Phil-P

Catalog #	Description	Dimension	Separation Range
SADP1002508 SADP1003008 SADP100508 SADP100308	AppliChrom ABOA DMSO-Phil-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SADP2002508 SADP2003008 SADP200508 SADP200308	AppliChrom ABOA DMSO-Phil-P-200	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-20.000Da
SADP2502508 SADP2503008 SADP250508 SADP250308	AppliChrom ABOA DMSO-Phil-P-250	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-70.000Da
SADP3002508 SADP3003008 SADP300508 SADP300308	AppliChrom ABOA DMSO-Phil-P-300	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-300kDa
SADP3502508 SADP3503008 SADP350508 SADP350308	AppliChrom ABOA DMSO-Phil-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-1.MioDa
SADP4002508 SADP4003008 SADP400508 SADP400308	AppliChrom ABOA DMSO-Phil-P-400	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10kDa-5 MioDa
SADP4502508 SADP4503008 SADP450508 SADP450308	AppliChrom ABOA DMSO-Phil-P-450	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100 kDa→10MioDa
SADPM2508 SADPM3008 SADPM508 SADPM308	AppliChrom ABOA DMSO-Phil-P-Multipore	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da - 1MioDa




Lignin

AppliChrom ABOA SuperOH-P

- optimized for aqueous GPC/SEC-separations
- aqueous high hydrophilic polymeric base material
- low column bleeding for low detector noise
- 7 μ particle technology (standard) for high platecounts and high resolution^{*)}
- large pore volume for high resolution
- pH stable 2.5-12
- pressure stability 30-80bar (depending on poresize)
- temperature stability 10-85°C
- individual poresizes for individual molecular weight ranges
- multi-pore technology for broad range of molecular size
- increased lifetime of GPC/SEC columns by combination of proprietary particle- and packing technology
- service-applicationcenter for methodscreening available in Oranienburg (Germany / Europe)

^{*)} 7 μ Particletechnology is standard for the small porous series – 100 and – 200.


GPC Calibration Curves AppliChrom ABOA SuperOH-P




Analyte: PEO/PEG

Column: AppliChrom ABOA SuperOH-P-**100**
 AppliChrom ABOA SuperOH-P-**200**
 AppliChrom ABOA SuperOH-P-**300**
 AppliChrom ABOA SuperOH-P-**400**

Dimension: e.a. 300mm x 8mm
 Mobil Phase: H₂O,
 Flow: 0,5ml/min
 Temperature: 20°C
 Detection: RI
 Injection: 20 μ l sample




For aqueous applications



Analyte: Oligosaccharide


Column: AppliChrom ABOASuperOH-P-100
Dimension: 2x 300mm x 8mm
Mobil Phase: H₂O
Flow: 0,5ml/min
Temperature: 40°C
Detection: RI
Injection: 20µl sample



Analyte: Pectin

sample with high content of oligomers


Column: AppliChrom ABOASuperOH-P-100
AppliChrom ABOASuperOH-P-350
Dimension: e.a. 300mm x 8mm
Mobil Phase: H₂O, NaKHPO₄ (pH 6.8 0.07M) + 50mM NaCl
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample



Analyte: PEGylated protein

Separation of an approx. 5.000Da product of a PEGylated protein of approx. 100.000dalton.


Column: AppliChrom ABOASuperOH-P-100
AppliChrom ABOASuperOH-P-350
Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,05% NaN₃ in H₂O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample



Analyte: Starch hydrolysate

analysis of a 100kDa fraction

Column: AppliChrom ABOASuperOH-P-100
AppliChrom ABOASuperOH-P-350
Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,2M NaN₃ in H₂O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample



Analyte:

Inulin

analysis of a 100Da- 1Mio Da,
including oligomer separation

Column:

AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

0,075M NaNO₃, 5g/l Na₂HPO₄·x7H₂O in H₂O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample



Analyte:

Hyaluronic acid

including oligomer separation

Column:

AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

0,075M NaNO₃, 5g/l Na₂HPO₄·x7H₂O in H₂O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample



Analyte:

Corn syrup

analysis of a 100Da- 1Mio Da,

Column:

AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

H₂O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample



Analyte:

**Water-based acrylate
dispersion**

analysis of a 100Da- 1Mio Da area

Column:

AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

1,0M NaNO₃ in H₂O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample




Analyte: Heparin-Na, 8-25kDa

including oligomer separation

Column: AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350

Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,075M NaNO₃, 5g/l Na₂HPO₄·x7H₂O in H₂O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample




Analyte: Dextran sulfate-Na

analysis of a 100Da- 1Mio Da,

Column: AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350

Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,075M NaNO₃, 5g/l Na₂HPO₄·x7H₂O in H₂O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample




Analyte: Alginate-Na

analysis of a 100Da- 1Mio Da,

Column: AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350

Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,075M NaNO₃, 5g/l Na₂HPO₄·x7H₂O in H₂O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample




Analyte: Carrageenan

analysis of a 1.000Da- 5Mio Da,

Column: AppliChrom ABOASuperOH-P-350
AppliChrom ABOASuperOH-P-450

Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,075M LiNO₃ in H₂O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample



Analyte:

Pullulan

including oligomer separation

Column:

AppliChrom ABOASuperOH-P-Screening

Dimension:

300mm x 8mm

Mobil Phase:

0,075M NaNO₃ in H₂O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample



Analyte:

Maltodextrin 12

analysis of a 100Da-1Mio Da,

Column:

AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

0,075M NaNO₃ in H₂O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample



Analyte:

Fish protein hydrolysates

3 different batches
analysis of a 100Da-70.000 Da,
including oligomer separation

Column:

AppliChrom ABOASuperOH-P-250

Dimension:

3x 300mm x 8mm

Mobil Phase:

0,05M NaNO₃ + 0,07M Na_{1,5}H_{1,5}PO₄ in H₂O

Flow:

1,0ml/min

Temperature:


30°C

Detection:

RI

Injection:

20µl sample



Analyte:

Polycarboxylate ether

3 different batches

Column:

AppliChrom ABOASuperOH-P-350

Dimension:

3x 300mm x 8mm

Mobil Phase:

0,1M NaNO₃ + 0,05M Na₂HPO₄ in H₂O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:


20µl sample



Analyte: Water-soluble casein fraction

Column: AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350


Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,01M $\text{Na}_{1,5}\text{H}_{1,5}\text{PO}_4$ & 0,3M NaCl in H_2O
Flow: 1,0ml/min
Temperature: 30°C
Detection: RI
Injection: 50µl sample



Analyte: Jelly Bean „Gummibärchen“

Column: AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350


Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,01M $\text{Na}_{1,5}\text{H}_{1,5}\text{PO}_4$ & 0,3M NaCl in H_2O
Flow: 1,0ml/min
Temperature: 50°C
Detection: RI
Injection: 20µl sample
Ingredient according to label: 6.9% protein (gelatin) dissolved in the eluent



**Analyte: Polyvinyl alcohols,
88% degree of hydrolysis**
3 different batches
including oligomer separation

Column: AppliChrom ABOASuperOH-P-350


Dimension: 2x 300mm x 8mm
Mobil Phase: 0,05M Na_2HPO_4 + 0,1M NaNO_3 in H_2O
Flow: 1,0ml/min
Temperature: 30°C
Detection: RI
Injection: 20µl sample



Analyte: PVP, polyvinyl pyrrolidone
3 different batches
including oligomer separation

Column: AppliChrom ABOASuperOH-P-350


Dimension: 2x 300mm x 8mm
Mobil Phase: 0,1M NaNO_3 in H_2O + 20% ACN
Flow: 1,0ml/min
Temperature: 30°C
Detection: RI
Injection: 50µl sample
PVP (M = 8,3 / 23,4 / 33,7 / 175kDa)
dissolved in the eluent



Analyte: **Pork gelatin vs. gelatin from collagen hydrolysate**
including oligomer separation

Column: AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350


Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,01M $\text{Na}_{1,5}\text{H}_{1,5}\text{PO}_4$ & 0,3M NaCl in H_2O
Flow: 1,0ml/min
Temperature: 30°C
Detection: RI
Injection: ea. 50µl sample
Pork gelatine (blue- - - -),
Gelatin from collagen hydrolyzate (red ----)



Analyte: **Pork gelatin**
analysis of a 100Da- 1Mio Da,

Column: AppliChrom ABOASuperOH-P-250
AppliChrom ABOASuperOH-P-350


Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,075M NaNO_3 , 5g/l $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$ in H_2O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample



GPC/SEC – Separation by size (ΔS) || **HPLC** – Separation according to bond strength (ΔH)
K = Partition coefficient

AppliChrom ABOA SuperOH-P

Catalog #	Description	Dimension	Separation Range
SASOH1002508 SASOH1003008 SASOH100508 SASOH100308	AppliChrom ABOA SuperOH-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SASOH2002508 SASOH2003008 SASOH200508 SASOH200308	AppliChrom ABOA SuperOH-P-200	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-20.000Da
SASOH2502508 SASOH2503008 SASOH250508 SASOH250308	AppliChrom ABOA SuperOH-P-250	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-70.000Da
SASOH3002508 SASOH3003008 SASOH300508 SASOH300308	AppliChrom ABOA SuperOH-P-300	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-300kDa
SASOH3502508 SASOH3503008 SASOH350508 SASOH350308	AppliChrom ABOA SuperOH-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	2.500Da-1.MioDa
SASOH4002508 SASOH4003008 SASOH400508 SASOH400308	AppliChrom ABOA SuperOH-P-400	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10kDa-5 MioDa
SASOH4502508 SASOH4503008 SASOH450508 SASOH450308	AppliChrom ABOA SuperOH-P-450	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	50 kDa→10MioDa
SASOHS2508 SASOHS3008 SASOHS508 SASOHS308	AppliChrom ABOA SuperOH-P-Screening	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da - 1MioDa



For HPLC-Analyses of


- Sugars/Carbohydrates
 - Sugar Alcohols
 - Alcohols
- Carboxylic Acids

Special Polymer for fast, easy and reliable determination using HPLC-RI or HPLC-ELSD at 60-80°C.


Advantages:

- Low cost for eluent supply. Eluent = water.
- Environmental friendly because: Eluent = water.
- Low cost for used eluent waste: Eluent = water.
- Easy to handle analysis. Direct analysis from aqueous sample!
- Low invest: Measurement possible with standard HPLC-system if connected with RI or ELSD (evaporative lightscattering detector).

- | | |
|---------------------------------------|---|
| AppliChrom ABOA SugarSep-Ca | - Analysis of sugars, sugar alcohols, alcohols. |
| AppliChrom ABOA SugarSep-Pb | - Analysis of sugars |
| AppliChrom ABOA SugarSep-H | - Analysis of sugars, sugaralcohols, alcohols and carboxylic acids. |
| AppliChrom ABOA SugarSep-Na | - Analysis of sugars, sugaralcohols, alcohols and carboxylic acids |
| AppliChrom ABOA SugarSep-Oligo | - Analysis of sugars, sugaralcohols, alcohols and carboxylic acids. |




Identification of many sugars

**Analyte:****Mixture of**

1. Sucrose (Saccarose)
2. Glucose
3. Fructose
4. Glycerin
5. Mannitol
6. Sorbitol

Column:
Dimension:
Mobil Phase:
Flow:
Temperature:
Detection:
Injection:


AppliChrom ABOA SugarSep-Ca
300mm x 8mm
H₂O
0,5ml/min
80°C
RI
20µl sample

**Analyte:****Mixture of**

1. Sucrose (Saccarose)
2. Glucose
3. Fructose
4. Glycerin

Column:
Dimension:
Mobil Phase:
Flow:
Temperature:
Detection:
Injection:


AppliChrom ABOA SugarSep-Pb
300mm x 8mm
H₂O
0,4ml/min
60°C
RI
20µl sample

**Analyte:****Mixture of**

1. Sucrose (Saccarose)
2. Glucose
3. Fructose
4. Glycerin
5. Ethanol

Column:
Dimension:
Mobil Phase:
Flow:
Temperature:
Detection:
Injection:


AppliChrom ABOA SugarSep-Ca
300mm x 8mm
H₂O
0,5ml/min
80°C
RI
20µl sample

**Analyte:****Honey (solvet in H₂O)**


1. Dp 3
2. Dp 2
3. Glucose
4. Fructose

Column:
Dimension:
Mobil Phase:
Flow:
Temperature:
Detection:
Injection:

AppliChrom ABOA SugarSep-Ca
300mm x 8mm
H₂O
0,5ml/min
80°C
RI
20µl sample

**Analyte: Honey**

Column: AppliChrom ABOA SugarSep-Na
 Dimension: 300mm x 8mm
 Mobil Phase: H₂O
 Flow: 0,5ml/min
 Temperature: 80°C
 Detection: RI
 Injection: 20µl sample


**Analyte: Malt Beer**

Column: AppliChrom ABOA SugarSep-Na
 Dimension: 300mm x 8mm
 Mobil Phase: H₂O
 Flow: 0,5ml/min
 Temperature: 80°C
 Detection: RI
 Injection: 20µl sample


AppliChrom BPT Technology: a combination of small, medium and large pores in each particle ensures an increase of calibration range. No matching porosities effects – respective artificial shoulders in chromatogrammes known from many column combinations are significant reduced – for improving your GPC Chromatography.

Conventional GPC technology: combining columns of various poresize or by combining different poresizes in one column enlarges the calibration range – but it can lead to artefacts in the exact calibration of the GPC system that reduces molecular size accuracy determination

Flowdirection (BPT technology)




Flowdirection (conventional mixed technology)




AppliChrom ABOA SugarSep

Catalog #	Description	Dimension	
SASCAI102508 SASCAI103008 SASCAI10508 SASCAI10308	AppliChrom ABOA SugarSep-Ca I	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASPBW102508 SASPBW103008 SASPBW10508 SASPBW10308	AppliChrom ABOA SugarSep-Pb-Wood	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASHI102508 SASHI103008 SASHI10508 SASHI10308	AppliChrom ABOA SugarSep-H I (SO ₃ H)	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASHII102508 SASHII103008 SASHII10508 SASHII10308	AppliChrom ABOA SugarSep-H II (SO ₃ H)	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASHIII102508 SASHIII103008 SASHIII10508 SASHIII10308	AppliChrom ABOA SugarSep-H III (SO ₃ H)	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASNA102508 SASNA103008 SASNA10508 SASNA10308	AppliChrom ABOA SugarSep-Na	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	




D-Glucose



L-Glucose


Oligosaccharide analysis in water



Analyte: Maltodextrin 19

[Detailed view](#)


Column: AppliChrom ABOA SugarSep-Oligo-Na
 Dimension: 300mm x 8mm
 Mobil Phase: H₂O
 Flow: 0,25ml/min
 Temperature: 70°C
 Detection: RI
 Injection: 20µl sample



Analyte: Maltodextrin 12

[Detailed view](#)


Column: AppliChrom ABOA SugarSep-Oligo-Na
 Dimension: 300mm x 8mm
 Mobil Phase: H₂O
 Flow: 0,25ml/min
 Temperature: 70°C
 Detection: RI
 Injection: 20µl sample



Analyte: Maltodextrin 6

[Detailed view](#)

Column: AppliChrom ABOA SugarSep-Oligo-Na
 Dimension: 300mm x 8mm
 Mobil Phase: H₂O
 Flow: 0,25ml/min
 Temperature: 70°C
 Detection: RI
 Injection: 20µl sample




Analyte: Maltodextrin 12 (black)

[Inulin \(blue\)](#)

Column: AppliChrom ABOA SugarSep-Oligo-Na
 Dimension: 300mm x 8mm
 Mobil Phase: H₂O
 Flow: 0,25ml/min
 Temperature: 70°C
 Detection: RI
 Injection: 20µl sample

Oligosaccharide analysis in water



Analyte: Maltodextrin 12
+ extra Glucose

Detailed view

Column: AppliChrom ABOA SugarSep-Oligo-Ag

Dimension: 300mm x 8mm


Mobil Phase: H₂O

Flow: 0,25ml/min

Temperature: 80°C

Detection: RI

Injection: 20µl sample



Analyte: Maltodextrin 12
without extra Glucose

Detailed view

Column: AppliChrom ABOA SugarSep-Oligo-Ag

Dimension: 300mm x 8mm

Mobil Phase: H₂O


Flow: 0,25ml/min

Temperature: 80°C

Detection: RI

Injection: 20µl sample

Catalog #	Description	Dimension	
SASOLNa101508 SASOLNa102508 SASOLNa103008 SASOLNa10508 SASOLNa10308	AppliChrom ABOA SugarSep-Oligo-Na	150mmx8mm 250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	Oligosaccharide analysis in water
SASOLAg101508 SASOLAg102508 SASOLAg103008 SASOLAg10508 SASOLAg10308	AppliChrom ABOA SugarSep-Oligo-Ag	150mmx8mm 250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	Oligosaccharide analysis in water



Maltodextrin

AppliChrom ABOA StyDiViBe-P

AppliChrom GPC-columns for GPC analysis of organic molecules using (THF, toluene, chloroform)**.

Oligomers and polymers including the new GPC media line from **AppliChrom BPT* synthesis technology** for large range, high resolution separations with increased accuracy in calibration

- spherical high porous styrene-divinylbenzene GPC-media
- large molecular weight range: 100 → 10.000.000Da
- optimized for 1ml/min flowrate when using 8mm ID-columns
- high pressure stability of 150-50bar, depending on porosity
- high capacity from AppliChrom 8mm ID GPC columns
(5% more capacity than 7,8mm columns, 14% more capacity than 7,5mm ID columns)
- plus extra high pore volume from AppliChrom GPC synthesis technology for an extra increasing of peak capacity and resolution
- proprietary AppliChrom GPC column packing procedure for accurate peak performance, low back pressures and an extension of column lifetime
- long column lifetime for reduction of costs even at high throughput screening applications
- high level of reproducibility
- high purity of AppliChrom GPC particles and columns for pure GPC mechanisms, low signal noise and reduction of „systempeaks“ after GPC run

GPC-examples (THF): Amylose acetat, amylose propionat,
butyl rubber, cellulose diacetat, cellulosenitrat, polybutadiene, polycarbonate,
polyisoprene, PMMA (polymethylmethacrylate), propylenglycol, polystyrene,
polymethylstyrene, natural rubber, PVC (polyvinylchloride), polyvinylacetate, epoxid
resins, polyisocyanate, polyols,
polyurethans, plant oils/triglycerids/diglycerids,....

GPC-examples (toluene): Silicones, polydimethylsiloxan


GPC for epoxid resins, oligomers, isocyanates, PMMA / polymethylmethacrylate,
polyethylmethacrylate, PS/polystyrene, vegetable oils /triglycerides/diglycerides,...,
polybutadiene, polyisoprene, silicon / siliconoil / polydimethylsiloxane (in toluene),
PEG / polyethylenglycol, polypropylenoxide, polyethylenglycol-polypropylen glycol-
copolymer, PVC/polyvinylchloride, PU / polyurethane, celluloseacetate,
diallylphthalate, dialkylphthalate, alkyd resin e.g...

* BPT: Broad Poredistribution Technology – easy calibration curve for large range of molecular sizes; lowered calibration artefacts compared to single pore combination columns.

** Eluent THF, toluene or chloroform must be dry. If other eluents are planned – please ask us in advance. Please avoid: Eluent containing e.g. water, alcohols, acetonitril – can be relevant e.g. when using a combined HPLC/GPC system and/or in case a degasser is part of the chromatography system, avoid drying of column.

AppliChrom ABOA StyDiViBe molecular weight range and optimum range of molecular weights

	Range of molar mass	Optimum resolution
35A	100-2.500Da	<400Da
100A	100-10.000Da	ca. 800Da
500A	100-30.000Da	1.000-3.000Da
1.000A	100-70.000Da	3.000-10.000Da
1.500A – BPT	100-120.000Da ^{*)}	4.000–15.000Da
10 ⁵ A - BPT	500-1.500.000Da ^{*)}	10.000-150.000Da
10 ⁶ A - BPT	1000-4.000.000Da ^{*)}	20.000-400.000Da
10 ⁷ A - BPT	1000->10.000.000Da ^{*)}	30.000-2.000.000Da

 Molecular sizes range of AppliChrom ABOA StyDiViBe GPC-Serie in detail:
THF GPC-Calibration curves
**Analyte:**

Column:

PMMA Standards
 AppliChrom ABOA StyDiViBe-P-**35A**
 AppliChrom ABOA StyDiViBe-P-**100A**
 AppliChrom ABOA StyDiViBe-P-**500A**
 AppliChrom ABOA StyDiViBe-P-**1000A**
 AppliChrom ABOA StyDiViBe-P-**1500A**
 AppliChrom ABOA StyDiViBe-P-**10⁵A**
 AppliChrom ABOA StyDiViBe-P-**10⁶A**
 AppliChrom ABOA StyDiViBe-P-**10⁷A**

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

THF

Flow:

1,0ml/min

Temperature:

20°C

Detection:


RI

Injection:

100µl sample


To cover a large range of molecular sizes GPC columns of suited porosities can be combined. This is the typical state of the art in many GPC laboratories. Also now it is useful if a special range of molecular sizes should be zoomed by GPC. But – in some single cases this also results in more or less obvious inhomogenities of calibration curves that itself makes accurate mathematics of calibration challenging. To increase accuracy and to simplify the calibration AppliChrom developed the BPT* synthesis technology for significant lowering the inhomogenities phenomena and to improve the results.

* BPT: Broad Pore distribution Technology – easy calibration curve for large range of molecular sizes; lowered calibration artefacts compared to single pore combination columns.




Analyte: Polystyrene (PS)
Mp = 578Da

Column: AppliChrom ABOA StyDiV-Be-P-35Å
Dimension: 300mm x 8mm
Mobil Phase: THF
Flow: 1,0ml/min
Temperature: 20°C
Detection: UV 265nm
Injection: 20µl sample




Analyte: Polymethylmethacrylate (PMMA) oligomer

Column: AppliChrom ABOA StyDiV-Be-P-100Å
Dimension: 2x 300mm x 8mm
Mobil Phase: THF
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample




Analyte: Grapeoil (Triglyceride GPC)

Column: AppliChrom ABOA StyDiV-Be-P-100Å
Dimension: 2x 300mm x 8mm
Mobil Phase: THF
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample




Analyte: Polyethylene glycol (PEG)
Mp = 200Da


Column: AppliChrom ABOA StyDiV-Be-P-100Å
Dimension: 2x 300mm x 8mm
Mobil Phase: THF
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample



Pore to small



Pore perfect



Pore to large

Analyte: Bisphenol-A-Epichlorhydrin resin (I)

Column: AppliChrom ABOA StyDiV-Be-P-100Å

Dimension: 2x 300mm x 8mm

Mobil Phase: THF

Flow: 1,0ml/min

Temperature: 20°C

Detection: RI

Injection: 20µl sample

High resolving GPC up to 10.000Da.

Analyte: Bisphenol-A-Epichlorhydrin resin (I)

Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT

Dimension: 2x 300mm x 8mm

Mobil Phase: THF

Flow: 0,5ml/min

Temperature: 20°C

Detection: RI

Injection: 20µl sample

Good resolution, calibration range respective linearity in the range of 100-120.000Da, obtained from the special AppliChrom BPT-technology. No „surprising“ porosity artefacts from mixing particles with pores of different size for covering the full range of molecular sizes. Great resolution even if 8µl RI measuring cell is used.

Analyte: Bisphenol-A-Epichlorhydrin resin (I)

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 300mm x 8mm

Mobil Phase: THF


Flow: 1,0ml/min

Temperature: 20°C

Detection: RI


Injection: 20µl sample

Good linearity from monomer up to 1.5Mio Dalton, ideal for porosity gap artefact reduced GPC screening of large ranges of molecular weights combined with maintaining of oligomer resolution




Analyte: Polyol hard foam

Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT
Dimension: 2x 300mm x 8mm
Mobil Phase: THF
Flow: 0,5ml/min
Temperature: 45°C
Detection: RI
Injection: 20µl sample



Analyte: Polyol soft foam


Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT
Dimension: 2x 300mm x 8mm
Mobil Phase: THF
Flow: 0,5ml/min
Temperature: 45°C
Detection: RI
Injection: 20µl sample



Analyte: Aliphatic polyester

(adipate polyester)
including fingerprint


Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT
Dimension: 2x 300mm x 8mm
Mobil Phase: THF
Flow: 0,5ml/min
Temperature: 45°C
Detection: RI
Injection: 20µl sample



Analyte: Polyether polyol

Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT
Dimension: 2x 300mm x 8mm
1x 50mmx 8mm
Mobil Phase: THF
Flow: 0,5ml/min
Temperature: 45°C
Detection: RI
Injection: 20µl sample

High resolution



Analyte: PVC

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm


Mobil Phase: THF

Flow: 1,0ml/min

Temperature: 45°C

Detection: RI

Injection: 20µl sample



Analyte: PMMA

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm


Mobil Phase: THF

Flow: 1,0ml/min

Temperature: 45°C

Detection: RI

Injection: 20µl sample



Analyte: PMMA enlarged

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm


Mobil Phase: THF

Flow: 1,0ml/min

Temperature: 45°C

Detection: RI

Injection: 20µl sample



Analyte: Silicone

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm

Mobil Phase: Toluene


Flow: 1,0ml/min

Temperature: 20°C


Detection: RI

Injection: 20µl sample

Broad calibration range, no porosity artefacts observed


**Analyte: Polystyrene**

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT
 Dimension: 2x 300mm x 8mm
 Mobil Phase: THF
 Flow: 1,0ml/min
 Temperature: 45°C
 Detection: RI
 Injection: 20µl sample


**Analyte: Polystyrene enlarged**

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT
 Dimension: 2x 300mm x 8mm
 Mobil Phase: THF
 Flow: 1,0ml/min
 Temperature: 45°C
 Detection: RI
 Injection: 20µl sample

Ve [ml] vs M PMMA [Da]



Ve vs. M (PMMA)

**GPC calibration curve**

Column: AppliChrom ABOA StyDiV-Be-P-35Å
 Dimension: 300mm x 8mm
 Mobil Phase: THF
 Flow: 1,0ml/min
 Temperature: 45°C
 Detection: RI
 Injection: 20µl sample

GPC calibration curve, large pore volume plus low exclusion limit for high oligomer resolution even with low backpressure.

GPC calibration curve

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT
 Dimension: 2x 300mm x 8mm
 Mobil Phase: THF
 Flow: 1,0ml/min
 Temperature: 20°C
 Detection: RI
 Injection: 20µl sample

AppliChrom ABOA StDiViBe 10E5A BPT-Technology: Large calibration range respective good to calculate calibration curve from monomer up to 1.5Mio Dalton, ideal technological step ahead if a combination of porosities results in „artificial shoulders“ of calibration curves and/or in artificial shoulders of broad distributed molecular weight polymers. An easy to use tool to improve quality of results for your analytes.

AppliChrom ABOA StyDiViBe-P

Catalog #	Description	Dimension	Separation Range
ASDVBP352508 ASDVBP353008 ASDVBP35508 ASDVBP35308	AppliChrom ABOA StyDiViBe-P-35Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
ASDVBP1002508 ASDVBP1003008 ASDVBP100508 ASDVBP100308	AppliChrom ABOA StyDiViBe-P-100Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-10.000Da
ASDVBP5002508 ASDVBP5003008 ASDVBP500508 ASDVBP500308	AppliChrom ABOA StyDiViBe-P-500Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-30.000Da
ASDVBP10002508 ASDVBP10003008 ASDVBP1000508 ASDVBP1000308	AppliChrom ABOA StyDiViBe-P-1000Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-70.000Da
ASDVBP15002508 ASDVBP15003008 ASDVBP1500508 ASDVBP1500308	AppliChrom ABOA StyDiViBe-P-1500Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.00Da-120.000Da
ASDVBP40002508 ASDVBP40003008 ASDVBP4000508 ASDVBP4000308	AppliChrom ABOA StyDiViBe-P-10 ⁴ Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-400.000Da
ASDVBP50002508 ASDVBP50003008 ASDVBP5000508 ASDVBP5000308	AppliChrom ABOA StyDiViBe-P-10 ⁵ Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	500Da-1.5MioDa
ASDVBP60002508 ASDVBP60003008 ASDVBP6000508 ASDVBP6000308	AppliChrom ABOA StyDiViBe-P-10 ⁶ Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10.000Da - 4MioDa
ASDVBP70002508 ASDVBP70003008 ASDVBP7000508 ASDVBP7000308	AppliChrom ABOA StyDiViBe-P-10 ⁷ Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100.000Da - 10MioDa




Runge picture „homemade“

AppliChrom ABOA CatPhil-P

Special GPC / SEC columns for the aqueous analysis of neutral, anionic and additionally cationic polymers.

For GPC / SEC analysis of

- Polycations, polyammonium compounds
- polyamines (chitosans), polyethylenimines,
- PEGylated polyethyleneimines
- polysaccharides
- polyanions (heparins, pectins, ...)




Analyte: Poly DADMAC

Polydiallyldimethylammonium chloride, polyquaternium-6, Mw = 100-200kDa, CAS [26062-79-3]

Column: AppliChrom ABOA CatPhil-P-100
AppliChrom ABOA CatPhil-P-350


Dimension: e.a. 300mm x 8mm
Mobil Phase: NaNO_3 + 0,2% formic acid in H_2O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 100µl sample



Analyte: High molecular weight chitosan, poliglusam, polyglucosamine, CAS [9012-76-4]

Column: AppliChrom ABOA CatPhil-P-400

Dimension: 3x 300mm x 8mm
Mobil Phase: NaNO_3 + 0,2% formic acid in H_2O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 100µl sample




Analyte: Chitosan sulfate

Column: AppliChrom ABOA CatPhil-P-100
AppliChrom ABOA CatPhil-P-350

Dimension: e.a. 300mm x 8mm
Mobil Phase: 0,05M Na_2HPO_4 + 0,1M NaNO_3 in H_2O
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 100µl sample

For aqueous polycations including

**Analyte: PEI (polyethyleneimine)**

linear, PEI 150, CAS [9002-98-6]

Mn = 4x10⁴, Mw = 4.2x10⁴,Mp = 3.9x10⁴Da;**PEI (polyethyleneimine)**


linear, PEI 25, CAS [9002-98-6]

Mn = 2.1x10³, Mw = 2.9x10³,Mp = 2.15x10³ Da

Column: AppliChrom ABOA CatPhil-P-100
AppliChrom ABOA CatPhil-P-350


Dimension: e.a. 300mm x 8mm
Mobil Phase: 0.1M NaCl + 0.2% TFA in H₂O
Flow: 1.0ml/min
Temperature: 20°C
Detection: RI
Injection: 100µl sample

AppliChrom ABOA CatPhil-P Calibration curve
Porosity P-100, P-200, P-350 and P400,
based on dextrans

**Analyte: Dextrans**

Column: AppliChrom ABOA CatPhil-P-100
AppliChrom ABOA CatPhil-P-200
AppliChrom ABOA CatPhil-P-350
AppliChrom ABOA CatPhil-P-400


Dimension: e.a. 300mm x 8mm
Mobil Phase: H₂O,
Flow: 1.0ml/min
Temperature: 20°C
Detection: RI
Injection: 100µl sample



Polyethylenimine (PEI)


AppliChrom ABOA CatPhil-P

Catalog #	Description	Dimension	Separation Range
SACP1002508 SACP1003008 SACP100508 SACP100308	AppliChrom ABOA CatPhil-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SACP2002508 SACP2003008 SACP200508 SACP200308	AppliChrom ABOA CatPhil-P-200	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-20.000Da
SACP3502508 SACP3503008 SACP350508 SACP350308	AppliChrom ABOA CatPhil-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	2.500Da-1.000.000Da
SACP4002508 SACP4003008 SACP400508 SACP400308	AppliChrom ABOA CatPhil-P-400	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10.000Da-5.000.000Da



AppliChrom® ABOA CatPhil-P-JLJ series - additionally enhanced hydrophilicity


Special GPC / SEC columns for the aqueous analysis of neutral, anionic and additionally cationic polymers - JLJ series with further increased hydrophilicity for extended application range - also in pure aqueous eluents (calibration with dextran / pullulan and new: additionally with PEO / PEG or p-2-vinylpyridine possible).



Analyte: **Poly(2-vinylpyridine)**
CAS [25014-15-7], Mw = 40.000Da

Column: AppliChrom ABOA CatPhil-P-350-JLJ


Dimension: 300mm x 8mm
Mobil Phase: 0.1M NaCl + 0,2% TFA in H₂O,
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample



Analyte: **Polyquaternium-33**
very high molecular weight, copolymer of trimethylaminoethyl acrylate salt and acrylamide, CAS [69418-26-4] (red) vs. high molecular weight polyacrylamide, CAS [9003-05-8] 15Mio Da (blue, dashed)

Column: AppliChrom ABOA CatPhil-P-100-JLJ
AppliChrom ABOA CatPhil-P-350-JLJ
AppliChrom ABOA CatPhil-P-500-JLJ


Dimension: e.a. 300mm x 8mm
Mobil Phase: 0.1M NaCl + 0,2% TFA in H₂O,
Flow: 1,0ml/min
Temperature: 40°C
Detection: RI
Injection: 100µl sample



Analyte: **Poly DADMAC**
Polydiallyldimethylammonium chloride, polyquaternium-6
Mw = 100-200kDa, CAS [26062-79-3]

Column: AppliChrom ABOA CatPhil-P-100-JLJ
AppliChrom ABOA CatPhil-P-350-JLJ
AppliChrom ABOA CatPhil-P-500-JLJ

Dimension: e.a. 300mm x 8mm
Mobil Phase: 0.1M NaCl + 0,2% TFA in H₂O,
Flow: 1,0ml/min
Temperature: 40°C
Detection: RI
Injection: 100µl sample




Analyte: Poly(vinylimidazole),
CAS [25232-42-2]

Column: AppliChrom ABOA CatPhil-P-350-JLJ

Dimension: 300mm x 8mm
Mobil Phase: 0.1M NaCl + 0,2% TFA in H₂O,
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 20µl sample

AppliChrom ABOA CatPhil-P-JLJ Calibration curve
 Porosity P-100, P-200, P-350 and P400,
 based on dextrans



Analyte: PEO
 Dextran
 Polyacrylamide

Column: AppliChrom ABOA CatPhil-P-100-JLJ
 AppliChrom ABOA CatPhil-P-350-JLJ
 AppliChrom ABOA CatPhil-P-500-JLJ

Dimension: e.a. 300mm x 8mm
Mobil Phase: H₂O,
Flow: 1,0ml/min
Temperature: 20°C
Detection: RI
Injection: 100µl sample

AppliChrom ABOA CatPhil-P-JLJ

Catalog #	Description	Dimension	Separation Range
SACP100JLJ2508 SACP100JLJ3008 SACP100JLJ508 SACP100JLJ308	AppliChrom ABOA CatPhil-P-100-JLJ	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SACP350JLJ2508 SACP350JLJ3008 SACP350JLJ508 SACP350JLJ308	AppliChrom ABOA CatPhil-P-350-JLJ	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	2.500Da-1.000.000Da
SACP500JLJ2508 SACP500JLJ3008 SACP500JLJ508 SACP500JLJ308	AppliChrom ABOA CatPhil-P-500-JLJ	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10.000Da- ~50.000.000Da

A general guarantee for the functionality with all analytes can not be given. In individual cases it is always advisable to discuss the issue with our product specialists. An application laboratory in Oranienburg is also available for your service. Please contact us.

Errors and changes excepted. Products for research and development. Not tested for clinical, diagnostic or food applications.

AppliChrom® ABOA HFIP-Phil-P

Special GPC / SEC columns for the GPC analysis of HFIP-soluble polymers.

For GPC / SEC analyzes from *)

Polyesters (polybutylene terephthalate / PBT / CAS 24968-12-5,
polyethylene terephthalate / PET / CAS 25038-59-9, polylactide
PLA / CAS 26100-51-6)

Polyamide 6 / PA6 / polycaprolactam / CAS 25038-54-4, polyamide 6-6 /
PA6-6 / polyhexamethylene adipamide / CAS131-17-2, polyamide 6-10 /
PA6-10 / poly(hexamethylene sebacamide) / PA 6-10, CAS 9011-52-3

Other (paraformaldehyde / polyoxymethylene / POM / polyacetal / CAS
30525-89-4 polyethylenimine / PEI / poly (iminoethylene / polyaziridine /
CAS 9002-98-6)

HFIP GPC / SEC calibration versus

- PMMA / polymethylmethacrylate / CAS 901-14-7


HFIP Accessories:

5mM CF₃COOK or 0.075M CF₃COONa are added to the eluent HFIP
for the suppression of electrostatic interactions and for artefact-free GPC.

Benefits AppliChrom ABOA HFIP-Phil-P GPC columns series over standard columns

- very large pore volume for high GPC resolution
- high GPC resolution for oligomers / condensates of 100-70,000Da
- high resolution GPC separation for the range 100-800.000Da
- large areas with high linearity calibratable areas without porosity artifacts
- low bleed
- standard columns also used for GPC-LS and GPC viscosity


*) GPC / SEC of polylactide (PLA) - see also AppliChrom ABOA StyDiViBe in CHCl₃, for polyethylenimine (PEI)
- see also: AppliChrom ABOA CatPhil-P in H₂O



Analyte: **Polymethylmethacrylat**
 (PMMA), CAS 9011-14-7, CAS131-17-2
 M = 901.000, 96.760, 32.500, 3196, 540Da

Column: AppliChrom ABOA HFIP-Phil-P-350


Dimension: 2x 300mm x 8mm
 Mobil Phase: HFIP, 5mM CF₃COONa
 Flow: 0,5ml/min
 Temperature: 40°C
 Detection: RI
 Injection: 100µl sample



Analyte: **Polyethylenterephthalat**
 (PET), CAS 25038-59-9

Column: AppliChrom ABOA HFIP-Phil-P-350


Dimension: 2x 300mm x 8mm
 Mobil Phase: HFIP, 5mM CF₃COONa
 Flow: 0,5ml/min
 Temperature: 40°C
 Detection: RI
 Injection: 100µl sample



Analyte: **Polyethylenterephthalat**
 (PET), CAS 25038-59-9

Column: AppliChrom ABOA HFIP-Phil-P-350

Dimension: 2x 300mm x 8mm
 Mobil Phase: HFIP, 5mM CF₃COONa
 Flow: 0,5ml/min
 Temperature: 40°C
 Detection: RI
 Injection: 100µl sample



Analyte: **Polyamide 6,6 (PA6-6),**
 CAS131-17-2

Column: AppliChrom ABOA HFIP-Phil-P-350

Dimension: 2x 300mm x 8mm
 Mobil Phase: HFIP, 5mM CF₃COOK
 Flow: 0,5ml/min
 Temperature: 40°C
 Detection: RI
 Injection: 100µl sample; 1g/l

Good for many HFIP soluble polymers


AppliChrom ABOA HFIP-Phil-P

Catalog #	Description	Dimension	Separation Range
SAHFIP1002508 SAHFIP1003008 SAHFIP100508 SAHFIP100308	AppliChrom ABOA HFIP-Phil-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SAHFIP3502508 SAHFIP3503008 SAHFIP350508 SAHFIP350308	AppliChrom ABOA HFIP-Phil-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	up to 1.000.000 Da
SAHFIP5002508 SAHFIP5003008 SAHFIP500508 SAHFIP00308	AppliChrom ABOA HFIP-Phil-P-500	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10.000Da- ~50.000.000Da

Principle of GPC/SEC

In detail:

Small polymers intrude to nearly all pores
 Medium sized polymers intrude only to some pores
 Very large polymers are totally excluded from the pores
 => separation of polymers according to molecular size



AppliChrom Aceton-AQ-Phil-P GPC-columns for GPC analysis of organic molecules using (80%acetone, 20%water)*.

Organosolv Lignin GPC analysis are now very easy and reliable to perform (Patent pending).

- optimized for GPC analysis of organosolv lignins
- easy and reliable to handle in acetone/water/traces formic acid // 80/20/1 // v/v/v
- compatible with evaporative detection (ELSD, MS)*
- preparative GPC fractionation of organosolv lignin without salt possible
- recycling GPC respective peak recycling GPC with enormous separation efficiency possible for isolation of individual organosolv lignin substances in semipreparative scale
- molecular weight calibration vs. PEO/PEG
- spherical high porous polymeric GPC-media with no silanol activity for pure GPC
- large molecular weight range: 100 - 1.000.000Da
- high pressure stability of 100bar
- high peak and separation capacity
- AppliChrom GPC columns – innovations and quality Made in Germany.


*) ask your ELSD respective your MS-detector supplier for compatibility or your individual detector resistance versus acetone!

AppliChrom Aceton-AQ-Phil GPC Series:

molecular weight range and optimum range of molecular weights

Poresize	Separation range	Optimum resolution range
250	100-70.000Da	100-10.000Da
350	100-1.000.000Da	1000-400.000Da

* Eluent 80/20 Aceton/waterv/v, molecular weight calibration vs. PEO/PEG



Analyte: Organosolv Lignin GPC

(range: 100-1.500.000Da)

Column: 2x AppliChrom Aceton-AQ-Phil-P-250,


Dimension: ea. 300mm x 8mm
 Mobil Phase: acetone/H₂O/formic acid

80/20/1 v/v/v

Flow: 1,0ml/min

Temperature: 25°C


Detection: RI



Analyte: **PEO/PEG GPC**
(8.000Da, 1.000Da, 106 Da)

Column: 2x AppliChrom Aceton-AQ-Phil-P-250

Dimension: ea. 300mm x 8mm
Mobil Phase: acetone/H₂O 80/20 v/v
Flow: 1,0ml/min
Temperature: 25°C
Detection: RI




Analyte: **Organosolv Lignin GPC**

Column: 3x AppliChrom Aceton-AQ-Phil-P-350


Dimension: ea. 300mm x 8mm
Mobil Phase: acetone/H₂O/formic acid 80/20 /1 v/v/v
Flow: 1,0ml/min
Temperature: 25°C
Detection: RI

Catalog #	Description	Dimension	Separation Range
SAAPP2503008 SAAPP250508 SAAPP250308	GPC-Column AppliChrom Aceton-Aq-Phil-P-250	300mm x 8mm 50mm x 8mm 30Mm x 8mm	100 - 70.000Da
SAAPP3503008 SAAPP350508 SAAPP350308	GPC-Column AppliChrom Aceton-Aq-Phil-P-350	300mm x 8mm 50mm x 8mm 30Mm x 8mm	100 - 1.000.000Da




Organosolv Lignin

- Up to 540.000 plates/separation in liquid chromatography for highest level purity products.
- Semipreparative and preparative scale
- Low solvent consumption and savings in workforce time
- Easy and efficient to handle including service hotline and suited Peak-Recycling-GPC system
- Large range of solvent polarity




AppliChrom Peak-Recycling-GPC Series

- efficiency
- effectiveness
- maximum purity
- time savings
- low solvent consumption




Classical manual reinjected semiprep GPC

Column lenght: 1x 600mm x 20mm
 Separation: ~ 30.000plates / separation
 Mobile phase: 1310min x 10ml/min
 => 13100ml = 13,1l
 necessary for separation



AppliChrom Peak-Recycling-GPC

Column lenght: 3x 600mm x 20mm
 Separation: ~ 540.000plates / separation
 Mobile phase: approx 50min x 10ml =
 approx 500ml
 necessary for separation



AppliChrom-Columns are easy to handle

You connect the columns with the suited Peak-Recycling-GPC system and laboratory operator let system work until needed purification result is monitored to be sufficient. In case of questions AppliChrom customers are assisted by the AppliChrom team members with more than 20 years hands on experience in liquid chromatography and GPC techniques (ask info@applichrom.de for assistance).

Solvent polarity

Polar

Non-polar

DMF/DMAc/NMP > 90/10 CHCl₃/MeOH > THF > Acetic acid ethylester > CHCl₃ > CH₂Cl₂ > Toluene^{*)}

^{*)} In case of using polar and viscous DMF/NMP/DMAc – reduction of flowrate, heating of column and in many cases some salt in solvent is necessary got get separation. Please ask for columns prepacked with final solvent that should be used.

Solvent polarity:

A big variation of solvent polarities can be used for AppliChrom Peak-Recycling-GPC columns (polystyrene based series)



Column Dimension here: 600mm x 20mm and 50mm x 20mm

Product	Exclusion Limit [Da]	Preparative Column performance ^{*)}			Precolumns ^{*)} For 20, 32 and 40mm ID columns: 50mm x 20mm
		20mm ID x	32mm ID x	40mm ID x	
AppliChrom Peak-Recycling- GPC Column					
35A	2500	> 42.000	> 42.000	> 42.000	> 35.000
100A	5.000	> 42.000	> 42.000	> 42.000	> 35.000
500A	20.000	> 42.000	> 42.000	> 42.000	> 35.000
1.000A	70.000	> 42.000	> 42.000	> 42.000	> 35.000

^{*)} Guaranteed performance in plates/meter, typical performance > 50.000/m. In recycling. More than 1 Mio plates /separation possible!

Column ID	Regular Flow rate [ml/min] ^{*)}	Max. Flow rate [ml/min] ^{*)}
20mm	10	15
32mm	20	30
40mm	30-40	40

^{*)} In case of using polar and viscous solvents like DMF/NMP/DMAc – reduction of flowrate, heating of column and in many cases some salt in solvent is necessary for a good separation. Please ask for columns prepacked with final solvent that should be used.

AppliChrom Peak-Recycling-GPC Columns

Catalog #	Description	Dimension	
CRC3560020 PCRC355020	Column AppliChrom Peak-Recycling-GPC, 35A Precolumn AppliChrom Peak-Recycling-GPC, 35A	600mmx20mm 50mm x 20mm	Other dimensions available.
CRC10060020 PCRC1005020	Column AppliChrom Peak-Recycling-GPC, 100A Precolumn AppliChrom Peak-Recycling-GPC, 100A	600mmx20mm 50mm x 20mm	Other dimensions available.
CRC50060020 PCRC5005020	Column AppliChrom Peak-Recycling-GPC, 500A Precolumn AppliChrom Peak-Recycling-GPC, 500A	600mmx20mm 50mm x 20mm	Other dimensions available.
CRC100060020 PCRC10005020	Column AppliChrom Peak-Recycling-GPC, 35A Precolumn AppliChrom Peak-Recycling-GPC, 35A	600mmx20mm 50mm x 20mm	Other dimensions available.

AppliChrom ABOA GPC/SEC Series:

AppliChrom ABOA GPC columns are designed for high performance GPC separations. These columns contain high crosslinked polymeric particles of high reproducibility, good selectivity and high stability. Before shipment, each AppliChrom ABOA GPC column is tested and supplied with its individual test chromatogram and specifications. For long live time of AppliChrom ABOA GPC columns and for long term satisfaction it is important for customers to take care about some points that will be given in the following text. Due to several different specific applications of customers, the text can not be absolute – therefore, if you have any question, please do not hesitate asking the AppliChrom service at sales@applichrom.de and inform your local distributor.

After receiving of column:

- Control if there is any obvious physical damage from shipping
- Check if the column you got is the column you ordered
- All columns are shipped in the solvent that is documented in the column test certificate
- Verify the column quality

Please take into account that performance of column also depends on the chromatographic system at customers site.

Following parameters can have an influence the verification of column performance:

- Loop size, degree of loop filling, concentration of sample
- Capillary diameters and length in chromatographic system
- Temperature of column
- Effective detector cell volume
- Data acquisition rate
- Fitting connections

AppliChrom ABOA GPC column series connectors:

- In case a set of columns of different porosities plus eventual precolumn is ordered, customer will be also supplied with the suited column connectors, of 1/16" OD. The fittings have the 10-32 type thread and ferrules that are pre-fixed to ensure tightening without voidvolume

Tubing and fittings:

- 8mm ID AppliChrom ABOA GPC column series are designed to be connected with 1/16" OD capillary via 10-32er threading with HPLC-System
- If connecting chromatographic system with columns, please take care about the depth of capillary coming from system - connected by ferrule and fitting to columns. If distance between capillary end and ferrule is too low you get extra void volumes and loss of peak performance, if distance between capillary and ferrule is too long, system can leak or column inlet can be destroyed
- **See also: AppliChrom HPLC / GPC column user guide**

Installation of column

- Rinse chromatographic system (without column) with degassed and filtered eluent (here: DMSO eluent) with a moderate flowrate, check if there is no leak in the system, there is no blockade in the system and please make sure that detector will not be blocked or getting a to high pressure! Maybe laboratory must be continuous good tempered (warm enough) to be sure DMSO is still fluid and viscosity is low enough! Please ensure this temperature also at weekends and/or holidays. Please document system pressure for the later used flowrate, at e.g. 0.4ml/min. If salts are added to the eluent, please take care to avoid precipitation of salt in system – as a result of this irreversible damages can occur
- Connect column respective columns with capillary from injector. Take care columns are connected in flow direction. The outlet of column (of the final column if there is a set of column) should be connected direct with waste, bypassing the detector system
- If using more than 1 column (series of columns); take the columns in the following order:
inlet => precolumn, maincolumn of smallest porosity, maincolumn of medium porosity, maincolumn of largest porosity => outlet
- Heat the columns at least to 50-60°C or to about 80°C for at least 1h
- Start pumping with 0.1ml/min, check if connections are tight, observe the outlet of columns to make sure a clear eluent flows out. Pump at least 5min / column at this flowrate, document the observed pressure.
- Increase stepwise to the final flowrate, e.g. 0.4ml/min, document the backpressure
- If backpressure is constant and eluent leaving the outlet pipe is clear, connect with detector (please make sure that connector is warm enough for ensuring fluidity of DMSO)
- After a total approx 5-10 column volumes baseline should be stable and measurements can start.
- If system is not needed for some time (over night, weekends,...), please lower flowrate e.g. to 0.05 or 0.1ml min, but do not change the column temperature. Every cooling and heating process of column means stress for column
- Please make sure that there is not a column pressure drop of more than 80bar per column or 30bar per Precolumn
- Avoid fast pressure increase or fast pressure lowering of columns
- If you inject sample, please take care sample is clear and free of particles (e.g. by filtration) and please take the lowest concentration that is necessary to get the results you need. Especially if using / analysing reactive samples, column lifetime can be increased significantly by using low concentrations of sample. Samples should be diluted with eluent to avoid eventual extra-peaks

Important:

If you have any question, please do not hesitate contacting the AppliChrom service, e.g. at sales@applichrom.de

We want to get you as satisfied customer with good products plus good service.

Guarantee:

AppliChrom guarantees the quality according to its test certificate at the date of shipment. Column warranty is voided if end-fittings are removed. Applications at customers site are not in the responsibility of supplier.

All information on these sides are put together with the biggest care, has been worked on and been checked. No liability or guarantee can be still assumed for content and objective correctness as well as actuality and completeness. All information without guarantee.


Installation:

- a) Rinse first your HPLC / GPC system (suction tube, pump, injector incl. injection loop as well as capillary connections up to the column entrance) with the mobile phase, e.g., with 1ml / min, note the originating counterpressure.
- b) Connect the capillary to the column entrance with the HPLC / GPC column. Provide for the fact that the capillary is optimally connected with the column entrance. Between the capillary end and the column should be
 - no dead volume originate (by too short capillary piece) – leads to line widening
 - the capillary end also not longer as the admission within the HPLC / GPC-Column – can lead to leakage in the column head.

Please follow,

the space of the capillary ends can vary between different HPLC / GPC suppliers. To make sure that the correct space of the capillary end are reached, the adaptation of the ferrule is recommended in respective HPLC / GPC hardware.

Charting of the influence of the capillary length



- c) **Pump** the eluent by the HPLC / GPC system consisting of pump, injector, capillary between column and injector, column, e.g., with 1ml / min with 4.6 mm of columns ID. Let go the first eluent 2-4ml then the column leave directly in the waste - under avoidance of the detectors (is valid for columns to max. 4.6 mm ID and you let go lengths from 50-250mm), with 8 mm of columns ID and lengths from 250-300mm please 15-20ml eluent directly in the waste – under avoidance of the detectors. If you control and document, besides, please pressure, rinse the column in this manner so long to the measured pressure is steady and the eluent clearly of the column flows out.
- d) **Connect** the column end from that of the eluent resigned with the detector system. Guarantee that yourselves, besides, no significantly raised system counterpressure compared with the before measured value from c) puts (tip to blockage of the detector system, e.g., from dried up eluent, tough or crystallised solvents, mature test, blocked, narrowed or broken off capillaries...) - blockage of the detector system can lead to the destruction with pressure excess – use please only accordingly trained, attentive staff!
- e) **Use** the fitting and ferrule connection between column and detector system in the manner that a) **no dead volume** (line widening) and b) no excess end of the inside-recumbent capillary originates (can lead to the destruction of the column exit). Ideally is the new touch-down of a connection from ferrule and fitting or, e.g., the use of PEEK-screw connections at the column exit.
- f) **Document** at regular intervals the system pressure.

Please, use for AppliChrom-precursor cartridges only the AppliChrom-precursor cartridge holder fitting for it!

All information on these sides are put together with the biggest care, has been worked on and been checked. No liability or guarantee can be still assumed for content and objective correctness as well as actuality and completeness. All information without guarantee.

General terms and conditions

The prices are free, in EURO, without VAT, this is calculated separately.

The listed prices are valid in the Federal Republic of Germany, deviations may result from changes in the manufacturers' prices.

Our offer is a business-to-business (B2B) offer and also addressed to selected traders and freelancers as well as universities and other public clients but not to private customers.

Delivery terms:

Federal Republic of Germany: Delivery is not free (EXW Oranienburg, INCOTERMS 2010). We charge a shipping and packaging fee depending on the size and weight of the goods or by arrangement.

(EU and third countries): Delivery is carried out DDU / EXW Oranienburg, INCOTERMS 2010. If no specification of the buyer is made, we determine the freight carrier and charge the buyer the shipping costs as well as all additionally incurred costs.

Payment:

From invoice date within 30 days net or by appointment.

Complaints:

Are possible within 8 business days after receipt of the goods.

In any case, they must be notified in writing.

The return of the goods is only permitted with our consent.

In case of justified defects, we have the choice between exchange of the goods, rework or reimbursement of the purchase price.

Complaints do not entitle the buyer to suspend his payment.

Any further damages claims of the buyer are excluded, as long as they are not caused by gross negligence or intentional. In case of gross negligence or intentional breach of our contractual or statutory obligations, our liability is limited to the foreseeable damage; we do not assume any liability for any defects.

Claims under the Product Liability Act are not affected by the above regulation.

If the price offer is obviously based on an incorrect price due to transmission errors or error, we are entitled to make a subsequent price correction.

If the customer rejects this correction or if no agreement can be reached between the parties to the contract, we are entitled to a right of withdrawal at any time, as far as legally permissible.

Warranty:

Corresponds to legal requirements.

However, the guarantee / warranty of the manufacturer applies to commercial products.

The guarantee for the properties of the columns is limited to the conditions of the original test chromatogram.

AppliChrom ensures the quality of the AppliChrom HPLC / GPC columns under the following conditions:

- If a column does not function properly, AppliChrom will re-install the column or replace it without incurring any costs for the customer.
- In the case of return of columns, AppliChrom must first issue a return authorization. Defective products must be accompanied by a written declaration of the defect. A redemption only takes place under the following conditions:
 - All columns must be checked upon receipt and AppliChrom must have all defects within 14 days after receipt of the column.
- Unless otherwise agreed, the warranty period is limited to a maximum of 90 days after delivery date.
- The guarantee for the properties of the columns is limited to the conditions of the original test chromatogram.

To void the warranty / guarantee:

- Removing or unscrewing the end fittings of the HPLC/GPC columns automatically voids column warranty.
- Physical damage to the column due to misuse or mishandling, including mechanical shocks.
- Damage to the stationary phase or the base material due to handling in the case of incorrect chemical conditions, unsuited solvent, temperatures or pressure conditions.
- Failure due to high backpressures caused by improper solvent and/or sample filtration practices causing particulate build-up or precipitation in the column or end-fitting.
- Incorrect selection of packing material made by customer for their particular use or incompatibility of equipment, etc.

Use restriction

Only for R & D as well as laboratory, not tested for pharmacy or medical diagnostics.

Buyer will not use AppliChrom's name in any publicity or advertising without AppliChrom's prior written approval.

Proprietary right of reservation exists until full payment of the goods.

Only place of performance and jurisdiction is Oranienburg/Germany.

We reserve the right to make technical changes.

Our company does not participate in a consignment procedure.

Application	Page	Application	Page
Acrylate water-based (range: 100-1.000.000Da)	14	PEI, polyethyleneimine	35
Alginate-Na (range: 100-1.000.000Da)	15	PEO / PEG	12, 43
Aliphatic polyester	30	PEO, Dextrans, Polyacryamide, comparison of	38
Bisphenol-A-Epichlorohydrin (I) (pore to large)	29	PET, Polyethylene terephthalate (in HFIP)	40
Bisphenol-A-Epichlorohydrin (I) (pore to small)	29	PMMA, Polymethylmethacrylate (in HFIP)	40
Caramel color GPC (range: 100-1.500.000Da)	10	PMMA, Polymethylmethacrylate (in THF)	31
Carrageenan (range: 100-1.000.000Da)	15	PMMA, Polymethylmethacrylate, oligomer	28
Casein fraction water-soluble	17	Poly DADMAC	34
Chitosan sulfate	34	Poly DADMAC (comparison 3 samples)	37
Chitosan, high molecular weight	34	Poly(2-vinylpyridine), Mw = 40.000Da	37
Comparison of 3 samples, Poly DADMAC	37	Poly(N-isopropylacrylamide)	7
Comparison of PEO, Dextrans, Polyacryamide	38	Poly(vinylimidazole)	38
Corn syrup (range: 100-1.000.000Da)	14	Poly[2-(butenyl)2-oxazoline-co-ethylenimine]	8
Dextran 650	6	Polyacryamide, comparison of PEO, dextrans	38
Dextran sulfate-Na (range: 100-1.000.000Da)	15	Polyamide 6,6, (PA6-6) (in HFIP)	40
Dextrans	35	Polybutyl methacrylate/	
Dextrans, comparison of PEO, Polyacryamide	38	Poly(dimethylamino-ethylmethacrylate) M=22kDa	9
DP2, mixture of...4	21	Polycarboxylate ether	16
DP3, mixture of...4	21	Polydiallyldimethylammonium chloride,	
Ethanol, mixture of...3	21	polyquaternium-6,	34
Fish protein hydrolysates (range: 100-70.000Da)	16	Polyester, aliphatic	30
Fructose, mixture of...1*	21	Polyether polyol	30
Fructose, mixture of...2*	21	Polyethylene glycol (PEG)	28
Fructose, mixture of...3*	21	Polyethylene terephthalate, PET (in HFIP)	40
Fructose, mixture of...4*	21	Polyethyleneimine, PEI	35
Glucose, mixture of...1*	21	Polyglucosamine, high molecular weight	34
Glucose, mixture of...2*	21	Polyglusam, high molecular weight	34
Glucose, mixture of...3*	21	Polymethylmethacrylate, PMMA (in HFIP)	40
Glucose, mixture of...4*	21	Polyol soft foam	30
Glycerin, mixture of...1*	21	Polyquaternium-33	37
Glycerin, mixture of...2*	21	Polysaccharide	6
Glycerin, mixture of...3*	21	Polystyrene (PS), Mp = 578Da	28
Grapeoil (Triglyceride GPC)	28	Polystyrene, PS (in THF)	32
Hard Foam, polyol	30	Polyvinyl alcohol M=22kDa	8
Heparin.Na, 8-25kDa	15	Polyvinyl alcohols, 88% degree of hydrolysis	17
Honey	22	Polyvinyl pyridine degraded	8
Humic acids / humates (range: 100-1.500.000Da)	10	Polyvinyl pyridine fractions	8
Hyaluronic acid	14	Polyvinylchloride, PVC (in THF)	31
Inulin (range: 100-1.000.000Da)	14	Pork gelatin (range: 100-1.000.000Da)	18
Inulin comparison to Maltodextrin 12	24	Pork gelatin vs. Gelatin from collagen hydrolysate	18
Jelly Bean „Gumminbärchen“	17	PS, Polystyrene (in THF)	32
Lignin conifer bark	6	PS, Polystyrene, Mp = 578Da	28
Lignin, organosolv GPC	42, 43	Pullulan	16
Malt Beer	22	PVC, polyvinylchloride (in THF)	31
Maltodextrin 6	24	PVP, Polyvinyl pyrrolidone	17
Maltodextrin 6 detailed view	24	Silicone (in Toluene)	31
Maltodextrin 12	24	Soft foam, polyol	30
Maltodextrin 12 detailed view	24	Sorbitol, mixture of...1	21
Maltodextrin 12 (range: 100-1.000.000Da)	16	Soy protein GPC	
Maltodextrin 12 with extra Glucose	24	(range:100-1.00.00Da (based on dextrans))	9
Maltodextrin 12 without extra Glucose	24	Spruce bark extract	6
Maltodextrin 19	24	Starch hydrolysates	13
Maltodextrin 19 detailed view	24	Sucrose (saccharose), mixture of...1	21
Mannitol, mixture of...1	21	Sucrose (saccharose), mixture of...2	21
Manuka honey protein GPC		Sucrose (saccharose), mixture of...3	21
(range 100-1.000.000Da (based on dextrans))	9	Triglyceride GPC, Grapeoil	28
Melamin (range: 100-1.500.000Da)	10	UF-resin	7
MUF-resin	7	Water-based acrylate (range: 100-1.000.000Da)	14
Oligosaccharide	13	Water-soluble casein fraction	17
Organosolv Lignin GPC	42, 43		
PA6-6, Polyamide 6,6 (in HFIP)	40		
Pea protein GPC			
(range: 100-1.000.000Da (based on dextrans))	9		
Pea starch	7		
Pectin	13		
PEG / PEO	12, 43		
PEG, Polyethylene glycol	28		
PEGylated protein	13		

*

Mixture 1	Mixture 2	Mixture 3	Mixture 4
1. Sucrose	1. Sucrose	1. Sucrose	1. Dp 3
2. Glucose	2. Glucose	2. Glucose	2. Dp 2
3. Fructose	3. Fructose	3. Fructose	3. Glucose
4. Glycerin	4. Glycerin	4. Glycerin	4. Fructose
5. Mannitol		5.Ethanol	
6. Sorbitol			

