


# Polymer-based GPC / SEC / HPLC Columns Made by AppliChrom



## Product Catalog 2020

Polymers  
Sugars  
Biopolymers  
Degradation Products

## Table of Contents

I.	About Us	1
II.	AppliChrom ABOA DMAC-Phil-P	3
III.	AppliChrom ABOA DMSO-Phil-P	5
IV.	AppliChrom ABOA SuperOH-P	12
V.	AppliChrom ABOA SugarSep	20
VI.	AppliChrom ABOA SugarSep-Oligo	24
VII.	AppliChrom ABOA StyDiViBe-P	26
VIII.	AppliChrom ABOA CatPhil-P	34
IX.	AppliChrom ABOA CatPhil-P-JLJ	37
X.	AppliChrom HFIP-Phil-P	39
XI.	AppliChrom ABOA Series GPC/SEC Column Guide	42
XII.	AppliChrom HPLC / GPC Column User Guide	44
XIII.	AppliChrom General Terms and Conditions	45
XIV.	AppliChrom Application List of Application / Product Catalog 2020	46



## About Us

AppliChrom team and customers grow together by supporting chromatographers in actual challenging and often difficult separations. You can count on AppliChrom with its own production site for chromatography media, columns and applications and more than 50 years of HPLC/GPC experience at the traditional historical place where 1850 chromatography was first invented by Professor Runge (Dipl.-Ing. (FH) Susanne & Dr. Christian Dauwe).

## 13 Reasons for AppliChrom

### 1) We take it personally

Especially the support of our customers. From the analytical question up to the chromatographic process we discuss and accompany you on the way to your chromatographic result.

### 2) Application and Technical Support from one Source

As your chromatography partner we organize and look at the whole chromatographic process – also under optimal use of the HPLC/GPC/SEC-devices already in use.

### 3) Increasing Quality

You are looking for the best way to get your optimum chromatographic result in HPLC, GPC or SEC. We do not only offer high quality and highly reproducible products for high resolution and lifetime, but also a suitable consultation to your individual chromatographic task.

### 4) Saving Time

Chromatographers who start developing a method have to finally produce a reliable method that give complete information about all components of interest from a sample. Maybe the method has to be compatible to different kinds of chromatography systems or environments later – taking these aspects into account from the beginning saves much time and cost for our customers. Our product specialists are glad to assist you from the beginning of a project to make you successful.

### 5) Experience

AppliChrom's experience of own production of chromatography media, columns and applications (HPLC, GPC and SEC) allows us to support you with the background knowledge of the chromatography media producer with staff having more than 50 years of relevant experience.

### 6) Individual

From the beginning of your request we integrate your individual wishes, ideas and needs into the support to your optimal chromatographic solution. Resolution, selectivity, compatibility to existing systems, time saving or all of this – please ask us. AppliChrom has own application laboratories, own manufacturing capacities and own development unions to support you with your individual needs.

### 7) Innovative

AppliChrom chromatography media are typically based on special customers' needs that were covered via the AppliChrom development process unit finishing with the process validation to high end products that are available now for your reliable product solutions. Thus many very unique and powerful solutions are offered now for HPLC, SEC and GPC by AppliChrom.

### 8) Reproducibility

AppliChrom works at the headquarter Oranienburg has its own quality control laboratories, the production of chromatography columns in small and large series but also units to produce chromatography media in small scale and for larger bulk amounts. The quality of each batch and each column is controlled to ensure the highest level of reproducibility, column for column, batch to batch and year after year.

### 9) OEM business

AppliChrom supplies well-chosen OEM customers with individual chromatography bulk or column lines. The production also contains customized lines for special needs. You do not find these products necessarily in our public catalogue. If you have a special need – please speak with us.

### 10) Chromatography – Origin from Oranienburg, innovations from Oranienburg

Maybe it is interesting for you when and where chromatography was invented. It was in 1850 in Oranienburg where professor Runge invented, published and commercialized chromatography products first in the world – it was an early form of paper chromatography – he called it: „Bilder die sich selber malen“ known also as Runge pictures. These are the roots of chromatography and the heritage of Oranienburg that we develop here at Oranienburg to innovative high end products to support you with your specific chromatographic solution.

### 11) Made in Europe

AppliChrom brand products are quality products made in Germany, coming from Oranienburg located in the Berlin region. Surrounded by several universities and institutes dealing with life and material science we get a continuous input of actual scientific challenges to get inspirations for needed new high end products. Anyway we are glad to hear what AppliChrom can do for you.

### 12) Carefree

Chromatographic separation should not be an adventure for the customer. We serve you from your request with the needed individual assistance in the beginning up to the after sales service/consultation.

### 13) Satisfaction

At the end AppliChrom's goal is to only have satisfied customers. This is our endeavour.

AppliChrom's core competence is the development and production of high quality liquid chromatography media in analytical and semipreparative scale. For this reason AppliChrom is using different reactor lines for producing several unique kinds of chromatography media. Especially for chromatography (HPLC, SEC, GPC) of biopolymers, oligomers, degradation products, fermentation products AppliChrom serves its customers with a unique product line.

AppliChrom branded products are quality products „Made in Germany“, these products are developed and produced in Germany and are submitted to a strict high class quality control.

The headquarter of AppliChrom is located in Oranienburg near Berlin. Oranienburg is the town where Friedlieb Ferdinand Runge invented the early form of paper chromatography in 1850.

Analytical chromatography is a process of separation, quantification and identification of chemical or biological compounds. Chromatography is used e.g. for research and analysis of food, drinks, pharmaceuticals, metabolomic processes, material sciences, environment and chemicals.

With our more than 5.000 products we supply chromatography columns (HPLC, SEC, GPC) for all current applications. Therefore, numerous institutes, universities, manufacturers, as well as our other customers appreciate the high quality, lifetime, reproducibility and reliability of AppliChrom products.

We extend our program steadily with new own chromatographic developments which fulfill the special demands of our customers to secure an economic workflow at the customers' site.

We are also glad about your inquiry!

## Important facts about AppliChrom

**Founded:** 2006, first market presence in 2009

**Company founder:** Dipl.-Ing. (FH) Susanne Dauwe (Chemical engineer)

**Research, development and chromatography applications:** Dr. Christian Dauwe (Chemist)

#### Products:

Silica- and polymeric based liquid chromatography media (HPLC, GPC, SEC)

HPLC-columns – C8, C18, HILIC & hybrid, ion exchanger and ion exclusion columns and media

Special easy to use HPLC columns for sugars, alcohols, acids (HILIC, ion exclusion, ion exchangers)

SEC/GPC columns for water soluble biopolymers

SEC/GPC columns for organic soluble polymers



Monument of Dr. Runge in front of the former Laboratory in Oranienburg

## AppliChrom ABOA DMAc-Phil-P

AppliChrom GPC columns for GPC analyzes organic (DMAc, DMF and NMP)

Spherical porous polymeric GPC chromatography phases


- wide range of pore and particle sizes
- optimal for 0.5ml / min flow rate with 8mm ID columns (40-50°C)
- Pressure stability 50-150bar, depending on the porosity
- high resolution due to high pore volume
- long service life
- high reproducibility
- high purity of the GPC column for good interaction-free GPC
- Please avoid: drying out
- Molar mass range: 100-1,000,000Da
- Examples: PMMA, PAN, cellulose, DMF and DMAc soluble polymers

AppliChrom ABOA DMAc-Phil-P molar mass range:


P-100	100-2.500Da
P-200	100-20.000Da
P-250	100-70.000Da
P-250JLD*)	100-100.000Da
P-300	1000-300.000Da
P-350	1.500-1.000.000Da

\*) JLD Pore: extended linear range

To cover a very wide range of molecular sizes, GPC columns of appropriate porosity can be combined.



### GPC Calibration Curves AppliChrom ABOA DMAc-Phil



Long linear calibration range

#### Analyte:

#### PMMA Standards

Column:

AppliChrom ABOA DMAcPhil-P-**100**  
 AppliChrom ABOA DMAcPhil-P-**200**  
 AppliChrom ABOA DMAcPhil-P-**250**  
 AppliChrom ABOA DMAcPhil-P-**250-JLD**  
 AppliChrom ABOA DMAcPhil-P-**300**  
 AppliChrom ABOA DMAcPhil-P-**350**

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

H<sub>2</sub>O

Flow:

0,5ml/min

Temperature:

45°C

Detection:


RI

Injection:

20µl sample

## AppliChrom ABOA DMAC-Phil-P

Catalog #	Description	Dimension	Separation Range
SADCP1002508 SADCP1003008 SADCP100508 SADCP100308	AppliChrom ABOA DMAC-Phil-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SADCP2002508 SADCP2003008 SADCP200508 SADCP200308	AppliChrom ABOA DMAC-Phil-P-200	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-20.000Da
SADCP2502508 SADCP2503008 SADCP250508 SADCP250308	AppliChrom ABOA DMAC-Phil-P-250	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-70.000Da
SADCP2502508J SADCP2503008J SADCP250508J SADCP250308J	AppliChrom ABOA DMAC-Phil-P-250	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-100.000Da
SADCP3002508 SADCP3003008 SADCP300508 SADCP300308	AppliChrom ABOA DMAC-Phil-P-300	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-300kDa
SADCP3502508 SADCP3503008 SADCP350508 SADCP350308	AppliChrom ABOA DMAC-Phil-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-1.MioDa
SADCP4002508 SADCP4003008 SADCP400508 SADCP400308	AppliChrom ABOA DMAC-Phil-P-400	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10kDa-5 MioDa
SADCP4502508 SADCP4503008 SADCP450508 SADCP450308	AppliChrom ABOA DMAC-Phil-P-450	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100 kDa→10MioDa



Reproducibility

## AppliChrom ABOA DMSO-Phil-P

Special GPC/SEC media for fast, accurate and robust GPC-analysis in DMSO


For GPC / SEC analysis in DMSO, examples:

- amylose, amylopectin, starch
- urea-formaldehyd resins (UF-resins)
- melamin-urea-formaldehyd resins (MUF-resins)
- lignins, humic substances, humic acids, coniferous wood bark essences
- polysaccharide, polysaccharid derivatives
- poly(N-isopropylacrylamid) PNIPA
- poly-vinylpyridin
- calibration: pullulan, dextran, polyvinylpyridin et al.


## Advantages

- optimized for DMSO-GPC applications
- interactionfree pure GPC/SEC
- easy, reliable and robust GPC/SEC-calibration by dextrans, pullulans et al.
- low column bleeding for low detectornoise for improved lightscattering or viscosity detection
- 12 $\mu$  particle technology for low backpressure
- large porevolume and optimized mass transfer for polymers giving optimized resolution
- low costs caused by long livetime of column – result of combination of optimized proprietary particle and packingtechnology.

Measuring range of selected AppliChrom ABOA DMSO-Phil-P GPC/SEC-columns, 300x8mm, poresizes 100, 200, 250, 350, 400




Different pore sizes available



**Analyte: Lignin conifer bark extracts**

Column: AppliChrom ABOA DMSO-Phil-P-250  
AppliChrom ABOA DMSO-Phil-P-350


Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,5ml/min  
Temperature: 80°C  
Detection: RI



**Analyte: Spruce bark extract**

Column: AppliChrom ABOA DMSO-Phil-P-250

Dimension: 2 x 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,5ml/min  
Temperature: 80°C  
Detection: RI




**Analyte: Polysaccharide**

(M ca. 70kDa)

Column: AppliChrom ABOA DMSO-Phil-P-250  
AppliChrom ABOA DMSO-Phil-P-350

Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,5ml/min  
Temperature: 80°C  
Detection: RI




**Analyte: Dextran 650**

Dextran from Leuconostoc spp.,  
M = 450.000-650.000Da + fructose

Column: AppliChrom ABOA DMSO-Phil-P-200  
AppliChrom ABOA DMSO-Phil-P-250  
AppliChrom ABOA DMSO-Phil-P-350


Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,5ml/min  
Temperature: 80°C  
Detection: RI



**Analyte: Pea starch**

Column: AppliChrom ABOA DMSO-Phil-P-200  
AppliChrom ABOA DMSO-Phil-P-250  
AppliChrom ABOA DMSO-Phil-P-350


Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,5ml/min  
Temperature: 80°C  
Detection: RI

**Analyte: MUF-resin**

completely DMSO-soluble melamin-urea formaldehyd resin (MUF-resin)

Column: AppliChrom ABOA DMSO-Phil-P-200  
AppliChrom ABOA DMSO-Phil-P-250  
AppliChrom ABOA DMSO-Phil-P-350


Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,5ml/min  
Temperature: 80°C  
Detection: RI  
Injection: 50µl sample

**Analyte: UF-resin**

DMSO-soluble urea formaldehyd resin

Column: AppliChrom ABOA DMSO-Phil-P-200  
AppliChrom ABOA DMSO-Phil-P-350

Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,5ml/min  
Temperature: 80°C  
Detection: RI


**Analyte: Poly(N-isopropylacrylamide)**

Further denominations: PNIPA, PNIPAAm, NIPA, PNIPAA or PNIPAm.  
CAS [25189-55-3], formula: (C<sub>6</sub>H<sub>11</sub>NO)<sub>n</sub>

3 different PNIPA fractions

Column: AppliChrom ABOA DMSO-Phil-P-300

Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,5ml/min  
Temperature: 80°C  
Detection: RI

**Analyte:****Polyvinylpyridine degraded**

Further denominations: CAS 25014-15-7,  $(\text{C}_7\text{H}_7\text{N})_n$   
low molecular weight (oligomeric)  
polyvinylpyridin fraction

Column: AppliChrom ABOA DMSO-Phil-P-250


Dimension: ea. 300mm x 8mm

Mobil Phase: DMSO / 0,075M  $\text{NaNO}_3$

Flow: 0,4ml/min

Temperature: 50°C

Detection: RI

**Analyte:****Polyvinylpyridine fractions**

Further denominations: CAS 25014-15-7,  $(\text{C}_7\text{H}_7\text{N})_n$   
75.7kDa, 20.9kDa, 3.2kDa

Superposition of 3 different polyvinylpyridin fractions

Column: AppliChrom ABOA DMSO-Phil-P-250


Dimension: ea. 300mm x 8mm

Mobil Phase: DMSO / 0,075M  $\text{NaNO}_3$

Flow: 0,4ml/min

Temperature: 50°C

Detection: RI

**Analyte:****Polyvinylalcohol M=22kDa**

Column: AppliChrom ABOA DMSO-Phil-P-300


Dimension: 2x. 300mm x 8mm

Mobil Phase: DMSO / 0,075M  $\text{NaNO}_3$

Flow: 0,4ml/min

Temperature: 50°C

Detection: RI

**Analyte:****Poly[2-(butenyl)2-oxazoline-co-ethylenimine]**

**M = 50.000Da**

Column: AppliChrom ABOA DMSO-Phil-P-300


Dimension: ea. 300mm x 8mm

Mobil Phase: DMSO / 0,075M  $\text{NaNO}_3$

Flow: 0,4ml/min


Temperature: 50°C

Detection: RI




**Analyte:** Polybutyl methacrylate/  
Poly(dimethylamino-  
ethylmethacrylate) M=22kDa

Column: AppliChrom ABOA DMSO-Phil-P-300  
Dimension: 2x. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,4ml/min  
Temperature: 50°C  
Detection: RI




**Analyte:** Pea protein GPC  
covering the calibration range of 100-  
1.000.000Da (based on dextrans)

Column: AppliChrom ABOA DMSO-Phil-P-Multipore  
Dimension: 3x. 300mm x 8mm  
Mobil Phase: DMSO  
Flow: 0,4ml/min  
Temperature: 55°C  
Detection: UV 280nm  
Calibration: vs. Dextran, pullulan or protein




**Analyte:** Soy protein GPC  
covering the calibration range of  
100-1.000.000Da (based on dextrans)

Column: AppliChrom ABOA DMSO-Phil-P-Multipore  
Dimension: 3x. 300mm x 8mm  
Mobil Phase: DMSO  
Flow: 0,4ml/min  
Temperature: 55°C  
Detection: UV 280nm  
Calibration: vs. Dextran, pullulan or protein



**Analyte:** Manuka honey protein GPC  
covering the calibration range of  
100-1.000.000Da (based on dextrans)

Column: AppliChrom ABOA DMSO-Phil-P-Multipore  
Dimension: 3x. 300mm x 8mm  
Mobil Phase: DMSO  
Flow: 0,4ml/min  
Temperature: 40°C  
Detection: UV 270nm  
Calibration: vs. Dextran, pullulan or protein




**Analyte: Caramel color GPC**

(range: 100Da-1.500.000Da)

SEC / GPC comparison of 2 Caramel  
Coulor

Column: AppliChrom ABOA DMSO-Phil-P-100  
AppliChrom ABOA DMSO-Phil-P-350

Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,3ml/min  
Temperature: 60°C  
Detection: UV 280nm




**Analyte: Humic acids / humates GPC**

(range: 100-1.500.000Da)

SEC / GPC comparison GPC / SEC  
comparison of 2 different humic acids /  
humates

Column: AppliChrom ABOA DMSO-Phil-P-Multipore

Dimension: 3x 300mm x 8mm  
Mobile Phase: DMSO  
Flow: 0,4ml/min  
Temperature: 70°C  
Detection: UV 280nm




**Analyte: Melamin GPC**

(range: 100-1.500.000Da)


Column: AppliChrom ABOA DMSO-Phil-P-100  
AppliChrom ABOA DMSO-Phil-P-350

Dimension: ea. 300mm x 8mm  
Mobil Phase: DMSO / 0,075M NaNO<sub>3</sub>  
Flow: 0,4ml/min  
Temperature: 60°C  
Detection: RI



## AppliChrom ABOA DMSO-Phil-P

Catalog #	Description	Dimension	Separation Range
SADP1002508 SADP1003008 SADP100508 SADP100308	AppliChrom ABOA DMSO-Phil-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SADP2002508 SADP2003008 SADP200508 SADP200308	AppliChrom ABOA DMSO-Phil-P-200	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-20.000Da
SADP2502508 SADP2503008 SADP250508 SADP250308	AppliChrom ABOA DMSO-Phil-P-250	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-70.000Da
SADP3002508 SADP3003008 SADP300508 SADP300308	AppliChrom ABOA DMSO-Phil-P-300	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-300kDa
SADP3502508 SADP3503008 SADP350508 SADP350308	AppliChrom ABOA DMSO-Phil-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-1.MioDa
SADP4002508 SADP4003008 SADP400508 SADP400308	AppliChrom ABOA DMSO-Phil-P-400	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10kDa-5 MioDa
SADP4502508 SADP4503008 SADP450508 SADP450308	AppliChrom ABOA DMSO-Phil-P-450	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100 kDa→10MioDa
SADPM2508 SADPM3008 SADPM508 SADPM308	AppliChrom ABOA DMSO-Phil-P-Multipore	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da - 1MioDa




Lignin

## AppliChrom ABOA SuperOH-P

- optimized for aqueous GPC/SEC-separations
- aqueous high hydrophilic polymeric base material
- low column bleeding for low detector noise
- 7 $\mu$  particle technology (standard) for high platecounts and high resolution\*)
- large pore volume for high resolution
- pH stable 2.5-12
- pressure stability 30-80bar (depending on poresize)
- temperature stability 10-85°C
- individual poresizes for individual molecular weight ranges
- multi-pore technology for broad range of molecular size
- increased lifetime of GPC/SEC columns by combination of proprietary particle- and packing technology
- service-applicationcenter for methodscreening available in Oranienburg (Germany / Europe)

\*) 7 $\mu$  Particletechnology is standard for the small porous series – 100 and – 200.


### GPC Calibration Curves AppliChrom ABOA SuperOH-P




#### Analyte: PEO/PEG

Column: AppliChrom ABOA SuperOH-P-**100**  
 AppliChrom ABOA SuperOH-P-**200**  
 AppliChrom ABOA SuperOH-P-**300**  
 AppliChrom ABOA SuperOH-P-**400**

Dimension: e.a. 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O,  
 Flow: 0,5ml/min  
 Temperature: 20°C  
 Detection: RI  
 Injection: 20 $\mu$ l sample




For aqueous applications



**Analyte: Oligosaccharide**


Column: AppliChrom ABOASuperOH-P-100  
Dimension: 2x 300mm x 8mm  
Mobil Phase: H<sub>2</sub>O  
Flow: 0,5ml/min  
Temperature: 40°C  
Detection: RI  
Injection: 20µl sample



**Analyte: Pectin**

sample with high content of oligomers


Column: AppliChrom ABOASuperOH-P-100  
AppliChrom ABOASuperOH-P-350  
Dimension: e.a. 300mm x 8mm  
Mobil Phase: H<sub>2</sub>O, NaKHPO<sub>4</sub> (pH 6.8 0.07M) + 50mM NaCl  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample



**Analyte: PEGylated protein**

Separation of an approx. 5.000Da product of a PEGylated protein of approx. 100.000dalton.


Column: AppliChrom ABOASuperOH-P-100  
AppliChrom ABOASuperOH-P-350  
Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,05% NaN<sub>3</sub> in H<sub>2</sub>O  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample



**Analyte: Starch hydrolysate**

analysis of a 100kDa fraction

Column: AppliChrom ABOASuperOH-P-100  
AppliChrom ABOASuperOH-P-350  
Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,2M NaNO<sub>3</sub> in H<sub>2</sub>O  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample



**Analyte:**

**Inulin**

analysis of a 100Da- 1Mio Da,  
including oligomer separation

Column:

AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

0,075M NaNO<sub>3</sub>, 5g/l Na<sub>2</sub>HPO<sub>4</sub>·x7H<sub>2</sub>O in H<sub>2</sub>O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample



**Analyte:**

**Hyaluronic acid**

including oligomer separation

Column:

AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

0,075M NaNO<sub>3</sub>, 5g/l Na<sub>2</sub>HPO<sub>4</sub>·x7H<sub>2</sub>O in H<sub>2</sub>O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample



**Analyte:**

**Corn syrup**

analysis of a 100Da- 1Mio Da,

Column:

AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

H<sub>2</sub>O

Flow:

1,0ml/min

Temperature:


20°C

Detection:

RI

Injection:

20µl sample



**Analyte:**

**Water-based acrylate  
dispersion**

analysis of a 100Da- 1Mio Da area

Column:

AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

1,0M NaNO<sub>3</sub> in H<sub>2</sub>O

Flow:

1,0ml/min

Temperature:

20°C


Detection:

RI

Injection:

20µl sample






**Analyte: Heparin-Na, 8-25kDa**

including oligomer separation

Column: AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,075M NaNO<sub>3</sub>, 5g/l Na<sub>2</sub>HPO<sub>4</sub>·x7H<sub>2</sub>O in H<sub>2</sub>O  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample




**Analyte: Dextran sulfate-Na**

analysis of a 100Da- 1Mio Da,

Column: AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,075M NaNO<sub>3</sub>, 5g/l Na<sub>2</sub>HPO<sub>4</sub>·x7H<sub>2</sub>O in H<sub>2</sub>O  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample




**Analyte: Alginate-Na**

analysis of a 100Da- 1Mio Da,

Column: AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,075M NaNO<sub>3</sub>, 5g/l Na<sub>2</sub>HPO<sub>4</sub>·x7H<sub>2</sub>O in H<sub>2</sub>O  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample




**Analyte: Carrageenan**

analysis of a 1.000Da- 5Mio Da,

Column: AppliChrom ABOASuperOH-P-350  
AppliChrom ABOASuperOH-P-450


Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,075M LiNO<sub>3</sub> in H<sub>2</sub>O  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample



**Analyte:** Pullulan  
including oligomer separation

**Column:** AppliChrom ABOASuperOH-P-Screening


**Dimension:** 300mm x 8mm  
**Mobil Phase:** 0,075M NaNO<sub>3</sub> in H<sub>2</sub>O  
**Flow:** 1,0ml/min  
**Temperature:** 20°C  
**Detection:** RI  
**Injection:** 20µl sample



**Analyte:** Maltodextrin 12  
analysis of a 100Da-1Mio Da,


**Column:** AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

**Dimension:** e.a. 300mm x 8mm  
**Mobil Phase:** 0,075M NaNO<sub>3</sub> in H<sub>2</sub>O  
**Flow:** 1,0ml/min  
**Temperature:** 20°C  
**Detection:** RI  
**Injection:** 20µl sample




**Analyte:** Fish protein hydrolysates  
3 different batches  
analysis of a 100Da-70.000 Da,  
including oligomer separation

**Column:** AppliChrom ABOASuperOH-P-250  
**Dimension:** 3x 300mm x 8mm  
**Mobil Phase:** 0,05M NaNO<sub>3</sub> + 0,07M Na<sub>1,5</sub>H<sub>1,5</sub>PO<sub>4</sub> in H<sub>2</sub>O  
**Flow:** 1,0ml/min  
**Temperature:** 30°C  
**Detection:** RI  
**Injection:** 20µl sample



**Analyte:** Polycarboxylate ether  
3 different batches


**Column:** AppliChrom ABOASuperOH-P-350  
**Dimension:** 3x 300mm x 8mm  
**Mobil Phase:** 0,1M NaNO<sub>3</sub> + 0,05M Na<sub>2</sub>HPO<sub>4</sub> in H<sub>2</sub>O  
**Flow:** 1,0ml/min  
**Temperature:** 20°C  
**Detection:** RI  
**Injection:** 20µl sample



**Analyte: Water-soluble casein fraction**

Column: AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350


Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,01M  $\text{Na}_{1,5}\text{H}_{1,5}\text{PO}_4$  & 0,3M NaCl in  $\text{H}_2\text{O}$   
Flow: 1,0ml/min  
Temperature: 30°C  
Detection: RI  
Injection: 50µl sample



**Analyte: Jelly Bean „Gummibärchen“**

Column: AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350


Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,01M  $\text{Na}_{1,5}\text{H}_{1,5}\text{PO}_4$  & 0,3M NaCl in  $\text{H}_2\text{O}$   
Flow: 1,0ml/min  
Temperature: 50°C  
Detection: RI  
Injection: 20µl sample  
Ingredient according to label: 6.9% protein (gelatin) dissolved in the eluent



**Analyte: Polyvinyl alcohols,  
88% degree of hydrolysis**  
3 different batches  
including oligomer separation

Column: AppliChrom ABOASuperOH-P-350


Dimension: 2x 300mm x 8mm  
Mobil Phase: 0,05M  $\text{Na}_2\text{HPO}_4$  + 0,1M  $\text{NaNO}_3$  in  $\text{H}_2\text{O}$   
Flow: 1,0ml/min  
Temperature: 30°C  
Detection: RI  
Injection: 20µl sample



**Analyte: PVP, polyvinyl pyrrolidone**  
3 different batches  
including oligomer separation

Column: AppliChrom ABOASuperOH-P-350

Dimension: 2x 300mm x 8mm  
Mobil Phase: 0,1M  $\text{NaNO}_3$  in  $\text{H}_2\text{O}$  + 20% ACN  
Flow: 1,0ml/min  
Temperature: 30°C  
Detection: RI  
Injection: 50µl sample  
PVP (M = 8,3 / 23,4 / 33,7 / 175kDa)  
dissolved in the eluent



**Analyte:** **Pork gelatin vs. gelatin from collagen hydrolysate**  
including oligomer separation

**Column:** AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

**Dimension:** e.a. 300mm x 8mm


**Mobil Phase:** 0,01M  $\text{Na}_{1,5}\text{H}_{1,5}\text{PO}_4$  & 0,3M NaCl in  $\text{H}_2\text{O}$

**Flow:** 1,0ml/min

**Temperature:** 30°C

**Detection:** RI

**Injection:** ea. 50µl sample  
Pork gelatine (blue- - - -),  
Gelatin from collagen hydrolyzate (red ----)



**Analyte:** **Pork gelatin**  
analysis of a 100Da- 1Mio Da,

**Column:** AppliChrom ABOASuperOH-P-250  
AppliChrom ABOASuperOH-P-350

**Dimension:** e.a. 300mm x 8mm


**Mobil Phase:** 0,075M  $\text{NaNO}_3$ , 5g/l  $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$  in  $\text{H}_2\text{O}$

**Flow:** 1,0ml/min

**Temperature:** 20°C

**Detection:** RI


**Injection:** 20µl sample



**GPC/SEC** – Separation by size ( $\Delta S$ ) || **HPLC** – Separation according to bond strength ( $\Delta H$ )  
**K** = Partition coefficient

## AppliChrom ABOA SuperOH-P

Catalog #	Description	Dimension	Separation Range
SASOH1002508 SASOH1003008 SASOH100508 SASOH100308	AppliChrom ABOA SuperOH-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SASOH2002508 SASOH2003008 SASOH200508 SASOH200308	AppliChrom ABOA SuperOH-P-200	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-20.000Da
SASOH2502508 SASOH2503008 SASOH250508 SASOH250308	AppliChrom ABOA SuperOH-P-250	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-70.000Da
SASOH3002508 SASOH3003008 SASOH300508 SASOH300308	AppliChrom ABOA SuperOH-P-300	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.000Da-300kDa
SASOH3502508 SASOH3503008 SASOH350508 SASOH350308	AppliChrom ABOA SuperOH-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	2.500Da-1.MioDa
SASOH4002508 SASOH4003008 SASOH400508 SASOH400308	AppliChrom ABOA SuperOH-P-400	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10kDa-5 MioDa
SASOH4502508 SASOH4503008 SASOH450508 SASOH450308	AppliChrom ABOA SuperOH-P-450	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	50 kDa→10MioDa
SASOHS2508 SASOHS3008 SASOHS508 SASOHS308	AppliChrom ABOA SuperOH-P-Screening	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da - 1MioDa



## For HPLC-Analyses of


- Sugars/Carbohydrates
  - Sugar Alcohols
  - Alcohols
- Carboxylic Acids

Special Polymer for fast, easy and reliable determination using HPLC-RI or HPLC-ELSD at 60-80°C.


### Advantages:

- Low cost for eluent supply. Eluent = water.
- Environmental friendly because: Eluent = water.
- Low cost for used eluent waste: Eluent = water.
- Easy to handle analysis. Direct analysis from aqueous sample!
- Low invest: Measurement possible with standard HPLC-system if connected with RI or ELSD (evaporative lightscattering detector).

- |                                       |   |
|---------------------------------------|---|
| <b>AppliChrom ABOA SugarSep-Ca</b>    | - Analysis of sugars, sugar alcohols, alcohols.                     |
| <b>AppliChrom ABOA SugarSep-Pb</b>    | - Analysis of sugars  |
| <b>AppliChrom ABOA SugarSep-H</b>     | - Analysis of sugars, sugaralcohols, alcohols and carobxylic acids. |
| <b>AppliChrom ABOA SugarSep-Na</b>    | - Analysis of sugars, sugaralcohols, alcohols and carobxylic acids  |
| <b>AppliChrom ABOA SugarSep-Oligo</b> | - Analysis of sugars, sugaralcohols, alcohols and carobxylic acids. |



Identification of many sugars




**Analyte:**

**Mixture of**

1. Sucrose (Saccarose)
2. Glucose
3. Fructose
4. Glycerin
5. Mannitol
6. Sorbitol

Column: AppliChrom ABOA SugarSep-Ca  
 Dimension: 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O  
 Flow: 0,5ml/min  
 Temperature: 80°C  
 Detection: RI  
 Injection: 20µl sample




**Analyte:**

**Mixture of**

1. Sucrose (Saccarose)
2. Glucose
3. Fructose
4. Glycerin

Column: AppliChrom ABOA SugarSep-Pb  
 Dimension: 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O  
 Flow: 0,4ml/min  
 Temperature: 60°C  
 Detection: RI  
 Injection: 20µl sample




**Analyte:**

**Mixture of**

1. Sucrose (Saccarose)
2. Glucose
3. Fructose
4. Glycerin
5. Ethanol

Column: AppliChrom ABOA SugarSep-Ca  
 Dimension: 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O  
 Flow: 0,5ml/min  
 Temperature: 80°C  
 Detection: RI  
 Injection: 20µl sample




**Analyte:**


**Honey (solvet in H<sub>2</sub>O)**

1. Dp 3
2. Dp 2
3. Glucose
4. Fructose

Column: AppliChrom ABOA SugarSep-Ca  
 Dimension: 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O  
 Flow: 0,5ml/min  
 Temperature: 80°C  
 Detection: RI  
 Injection: 20µl sample

**Analyte: Honey**

Column: AppliChrom ABOA SugarSep-Na  
 Dimension: 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O  
 Flow: 0,5ml/min  
 Temperature: 80°C  
 Detection: RI  
 Injection: 20µl sample


**Analyte: Malt Beer**

Column: AppliChrom ABOA SugarSep-Na  
 Dimension: 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O  
 Flow: 0,5ml/min  
 Temperature: 80°C  
 Detection: RI  
 Injection: 20µl sample


AppliChrom BPT Technology: a combination of small, medium and large pores in each particle ensures an increase of calibration range. No matching porosities effects – respective artificial shoulders in chromatogrammes known from many column combinations are significant reduced – for improving your GPC Chromatography.

Conventional GPC technology: combining columns of various poresize or by combining different poresizes in one column enlarges the calibration range – but it can lead to artefacts in the exact calibration of the GPC system that reduces molecular size accuracy determination

Flowdirection (BPT technology)




Flowdirection (conventional mixed technology)






## AppliChrom ABOA SugarSep

Catalog #	Description	Dimension	
SASCAI102508 SASCAI103008 SASCAI10508 SASCAI10308	AppliChrom ABOA SugarSep-Ca I	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASPBW102508 SASPBW103008 SASPBW10508 SASPBW10308	AppliChrom ABOA SugarSep-Pb-Wood	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASHI102508 SASHI103008 SASHI10508 SASHI10308	AppliChrom ABOA SugarSep-H I (SO <sub>3</sub> H)	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASHII102508 SASHII103008 SASHII10508 SASHII10308	AppliChrom ABOA SugarSep-H II (SO <sub>3</sub> H)	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASHIII102508 SASHIII103008 SASHIII10508 SASHIII10308	AppliChrom ABOA SugarSep-H III (SO <sub>3</sub> H)	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	
SASNA102508 SASNA103008 SASNA10508 SASNA10308	AppliChrom ABOA SugarSep-Na	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	




D-Glucose




L-Glucose


## Oligosaccharide analysis in water

**Analyte: Maltodextrin 19**


Column: AppliChrom ABOA SugarSep-Oligo  
Dimension: 300mm x 8mm  
Mobil Phase: H<sub>2</sub>O  
Flow: 0,25ml/min  
Temperature: 70°C  
Detection: RI  
Injection: 20µl sample

**Analyte: Maltodextrin 19**  
**Detailed view**

Column: AppliChrom ABOA SugarSep-Oligo  
Dimension: 300mm x 8mm  
Mobil Phase: H<sub>2</sub>O  
Flow: 0,25ml/min  
Temperature: 70°C  
Detection: RI  
Injection: 20µl sample


**Analyte: Maltodextrin 12**

Column: AppliChrom ABOA SugarSep-Oligo  
Dimension: 300mm x 8mm  
Mobil Phase: H<sub>2</sub>O  
Flow: 0,25ml/min  
Temperature: 70°C  
Detection: RI  
Injection: 20µl sample


**Analyte: Maltodextrin 12**  
**Detailed view**

Column: AppliChrom ABOA SugarSep-Oligo  
Dimension: 300mm x 8mm  
Mobil Phase: H<sub>2</sub>O  
Flow: 0,25ml/min  
Temperature: 70°C  
Detection: RI  
Injection: 20µl sample

## Oligosaccharide analysis in water


**Analyte: Maltodextrin 6**

Column: AppliChrom ABOA SugarSep-Oligo  
 Dimension: 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O  
 Flow: 0,25ml/min  
 Temperature: 70°C  
 Detection: RI  
 Injection: 20µl sample



**Analyte: Maltodextrin 6**

Detail view

Column: AppliChrom ABOA SugarSep-Oligo  
 Dimension: 300mm x 8mm  
 Mobil Phase: H<sub>2</sub>O  
 Flow: 0,25ml/min  
 Temperature: 70°C  
 Detection: RI  
 Injection: 20µl sample

### Oligo SEC/GPC of dextran and inulin

Catalog #	Description	Dimension	
SASOL101508 SASOL102508 SASOL103008 SASOL10508 SASOL10308	AppliChrom ABOA SugarSep-Oligo	150mmx8mm 250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	Oligosaccharide analysis in water



Maltodextrin

## AppliChrom ABOA StyDiViBe-P

AppliChrom GPC-columns for GPC analysis of organic molecules using (THF, toluene, chloroform)\*\*.

Oligomers and polymers including the new GPC media line from **AppliChrom BPT\* synthesis technology** for large range, high resolution separations with increased accuracy in calibration

- spherical high porous styrene-divinylbenzene GPC-media
- large molecular weight range: 100 → 10.000.000Da
- optimized for 1ml/min flowrate when using 8mm ID-columns
- high pressure stability of 150-50bar, depending on porosity
- high capacity from AppliChrom 8mm ID GPC columns  
(5% more capacity than 7,8mm columns, 14% more capacity than 7,5mm ID columns)
- plus extra high pore volume from AppliChrom GPC synthesis technology for an extra increasing of peak capacity and resolution
- proprietary AppliChrom GPC column packing procedure for accurate peak performance, low back pressures and an extension of column lifetime
- long column lifetime for reduction of costs even at high throughput screening applications
- high level of reproducibility
- high purity of AppliChrom GPC particles and columns for pure GPC mechanisms, low signal noise and reduction of „systempeaks“ after GPC run

GPC-examples (THF): Amylose acetat, amylose propionat,  
butyl rubber, cellulose diacetat, cellulosenitrat, polybutadiene, polycarbonate,  
polyisoprene, PMMA (polymethylmethacrylate), propylenglycol, polystyrene,  
polymethylstyrene, natural rubber, PVC (polyvinylchloride), polyvinylacetate, epoxid  
resins, polyisocyanate, polyols,  
polyurethans, plant oils/triglycerids/diglycerids,....

GPC-examples (toluene): Silicones, polydimethylsiloxan


GPC for epoxid resins, oligomers, isocyanates, PMMA / polymethylmethacrylate,  
polyethylmethacrylate, PS/polystyrene, vegetable oils /triglycerides/diglycerides,...,  
polybutadiene, polyisoprene, silicon / siliconoil / polydimethylsiloxane (in toluene),  
PEG / polyethylenglycol, polypropylenoxide, polyethylenglycol-polypropylen glycol-  
copolymer, PVC/polyvinylchloride, PU / polyurethane, celluloseacetate,  
diallylphthalate, dialkylphthalate, alkyd resin e.g...

\* BPT: Broad Poredistribution Technology – easy calibration curve for large range of molecular sizes; lowered calibration artefacts compared to single pore combination columns.

\*\* Eluent THF, toluene or chloroform must be dry. If other eluents are planned – please ask us in advance. Please avoid: Eluent containing e.g. water, alcohols, acetonitril – can be relevant e.g. when using a combined HPLC/GPC system and/or in case a degasser is part of the chromatography system, avoid drying of column.

**AppliChrom ABOA StyDiViBe molecular weight range and optimum range of molecular weights**

	Range of molar mass	Optimum resolution
35A	100-2.500Da	<400Da
100A	100-10.000Da	ca. 800Da
500A	100-30.000Da	1.000-3.000Da
1.000A	100-70.000Da	3.000-10.000Da
1.500A – BPT	100-120.000Da <sup>*)</sup>	4.000–15.000Da
10 <sup>5</sup> A - BPT	500-1.500.000Da <sup>*)</sup>	10.000-150.000Da
10 <sup>6</sup> A - BPT	1000-4.000.000Da <sup>*)</sup>	20.000-400.000Da
10 <sup>7</sup> A - BPT	1000->10.000.000Da <sup>*)</sup>	30.000-2.000.000Da

**Molecular sizes range of AppliChrom ABOA StyDiViBe GPC-Serie in detail:  
THF GPC-Calibration curves**
**Analyte:**

Column:

**PMMA Standards**

AppliChrom ABOA StyDiViBe-P-**35A**  
 AppliChrom ABOA StyDiViBe-P-**100A**  
 AppliChrom ABOA StyDiViBe-P-**500A**  
 AppliChrom ABOA StyDiViBe-P-**1000A**  
 AppliChrom ABOA StyDiViBe-P-**1500A**  
 AppliChrom ABOA StyDiViBe-P-**10<sup>5</sup>A**  
 AppliChrom ABOA StyDiViBe-P-**10<sup>6</sup>A**  
 AppliChrom ABOA StyDiViBe-P-**10<sup>7</sup>A**

Dimension:

e.a. 300mm x 8mm

Mobil Phase:

THF

Flow:

1,0ml/min

Temperature:

20°C

Detection:


RI

Injection:

100µl sample


To cover a large range of molecular sizes GPC columns of suited porosities can be combined. This is the typical state of the art in many GPC laboratories. Also now it is useful if a special range of molecular sizes should be zoomed by GPC. But – in some single cases this also results in more or less obvious inhomogenities of calibration curves that itself makes accurate mathematics of calibration challenging. To increase accuracy and to simplify the calibration AppliChrom developed the BPT\* synthesis technology for significant lowering the inhomogenities phenomena and to improve the results.

\* BPT: Broad Pore distribution Technology – easy calibration curve for large range of molecular sizes; lowered calibration artefacts compared to single pore combination columns.




**Analyte:** Polystyrene (PS)  
**Mp = 578Da**

Column: AppliChrom ABOA StyDiV-Be-P-35Å  
Dimension: 300mm x 8mm  
Mobil Phase: THF  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: UV 265nm  
Injection: 20µl sample




**Analyte:** Polymethylmethacrylate (PMMA) oligomer

Column: AppliChrom ABOA StyDiV-Be-P-100Å  
Dimension: 2x 300mm x 8mm  
Mobil Phase: THF  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample




**Analyte:** Grapeoil (Triglyceride GPC)

Column: AppliChrom ABOA StyDiV-Be-P-100Å  
Dimension: 2x 300mm x 8mm  
Mobil Phase: THF  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample




**Analyte:** Polyethylene glycol (PEG)  
**Mp = 200Da**


Column: AppliChrom ABOA StyDiV-Be-P-100Å  
Dimension: 2x 300mm x 8mm  
Mobil Phase: THF  
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 20µl sample



Pore to small



Pore perfect



Pore to large

**Analyte: Bisphenol-A-Epichlorhydrin resin (I)**

Column: AppliChrom ABOA StyDiV-Be-P-100Å

Dimension: 2x 300mm x 8mm

Mobil Phase: THF

Flow: 1,0ml/min

Temperature: 20°C

Detection: RI

Injection: 20µl sample

High resolving GPC up to 10.000Da.

**Analyte: Bisphenol-A-Epichlorhydrin resin (I)**

Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT

Dimension: 2x 300mm x 8mm

Mobil Phase: THF

Flow: 0,5ml/min

Temperature: 20°C

Detection: RI

Injection: 20µl sample

Good resolution, calibration range respective linearity in the range of 100-120.000Da, obtained from the special AppliChrom BPT-technology. No „surprising“ porosity artefacts from mixing particles with pores of different size for covering the full range of molecular sizes. Great resolution even if 8µl RI measuring cell is used.

**Analyte: Bisphenol-A-Epichlorhydrin resin (I)**

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 300mm x 8mm

Mobil Phase: THF


Flow: 1,0ml/min

Temperature: 20°C

Detection: RI


Injection: 20µl sample

Good linearity from monomer up to 1.5Mio Dalton, ideal for porosity gap artefact reduced GPC screening of large ranges of molecular weights combined with maintaining of oligomer resolution




**Analyte: Polyol hard foam**

Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT  
Dimension: 2x 300mm x 8mm  
Mobil Phase: THF  
Flow: 0,5ml/min  
Temperature: 45°C  
Detection: RI  
Injection: 20µl sample



**Analyte: Polyol soft foam**


Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT  
Dimension: 2x 300mm x 8mm  
Mobil Phase: THF  
Flow: 0,5ml/min  
Temperature: 45°C  
Detection: RI  
Injection: 20µl sample



**Analyte: Aliphatic polyester**

(adipate polyester)  
including fingerprint

Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT  
Dimension: 2x 300mm x 8mm  
Mobil Phase: THF  
Flow: 0,5ml/min  
Temperature: 45°C  
Detection: RI  
Injection: 20µl sample




**Analyte: Polyether polyol**

Column: AppliChrom ABOA StyDiV-Be-P-1500Å-BPT  
Dimension: 2x 300mm x 8mm  
1x 50mmx 8mm  
Mobil Phase: THF  
Flow: 0,5ml/min  
Temperature: 45°C  
Detection: RI  
Injection: 20µl sample

High resolution





**Analyte:** PVC

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm


Mobil Phase: THF

Flow: 1,0ml/min

Temperature: 45°C

Detection: RI

Injection: 20µl sample



**Analyte:** PMMA

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm


Mobil Phase: THF

Flow: 1,0ml/min

Temperature: 45°C

Detection: RI

Injection: 20µl sample



**Analyte:** PMMA enlarged

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm


Mobil Phase: THF

Flow: 1,0ml/min

Temperature: 45°C

Detection: RI

Injection: 20µl sample



**Analyte:** Silicone

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm

Mobil Phase: Toluene


Flow: 1,0ml/min

Temperature: 20°C

Detection: RI


Injection: 20µl sample

Broad calibration range, no porosity artefacts observed

**Analyte: Polystyrene**

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT


Dimension: 2x 300mm x 8mm  
 Mobil Phase: THF  
 Flow: 1,0ml/min  
 Temperature: 45°C  
 Detection: RI  
 Injection: 20µl sample

**Analyte: Polystyrene enlarged**


Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm  
 Mobil Phase: THF  
 Flow: 1,0ml/min  
 Temperature: 45°C  
 Detection: RI  
 Injection: 20µl sample

Ve [ml] vs M PMMA [Da]



Ve vs. M (PMMA)

**GPC calibration curve**

Column: AppliChrom ABOA StyDiV-Be-P-35Å

Dimension: 300mm x 8mm  
 Mobil Phase: THF  
 Flow: 1,0ml/min  
 Temperature: 45°C  
 Detection: RI  
 Injection: 20µl sample

GPC calibration curve, large pore volume plus low exclusion limit for high oligomer resolution even with low backpressure.

**GPC calibration curve**

Column: AppliChrom ABOA StyDiV-Be-P-10E5Å-BPT

Dimension: 2x 300mm x 8mm  
 Mobil Phase: THF  
 Flow: 1,0ml/min  
 Temperature: 20°C  
 Detection: RI  
 Injection: 20µl sample

AppliChrom ABOA StDiViBe 10E5A BPT-Technology: Large calibration range respective good to calculate calibration curve from monomer up to 1.5Mio Dalton, ideal technological step ahead if a combination of porosities results in „artificial shoulders“ of calibration curves and/or in artificial shoulders of broad distributed molecular weight polymers. An easy to use tool to improve quality of results for your analytes.

## AppliChrom ABOA StyDiViBe-P

Catalog #	Description	Dimension	Separation Range
ASDVBP352508 ASDVBP353008 ASDVBP35508 ASDVBP35308	AppliChrom ABOA StyDiViBe-P-35Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
ASDVBP1002508 ASDVBP1003008 ASDVBP100508 ASDVBP100308	AppliChrom ABOA StyDiViBe-P-100Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-10.000Da
ASDVBP5002508 ASDVBP5003008 ASDVBP500508 ASDVBP500308	AppliChrom ABOA StyDiViBe-P-500Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-30.000Da
ASDVBP10002508 ASDVBP10003008 ASDVBP1000508 ASDVBP1000308	AppliChrom ABOA StyDiViBe-P-1000Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-70.000Da
ASDVBP15002508 ASDVBP15003008 ASDVBP1500508 ASDVBP1500308	AppliChrom ABOA StyDiViBe-P-1500Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	1.00Da-120.000Da
ASDVBP4002508 ASDVBP4003008 ASDVBP400508 ASDVBP400308	AppliChrom ABOA StyDiViBe-P-10 <sup>4</sup> Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-400.000Da
ASDVBP50002508 ASDVBP50003008 ASDVBP5000508 ASDVBP5000308	AppliChrom ABOA StyDiViBe-P-10 <sup>5</sup> Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	500Da-1.5MioDa
ASDVBP60002508 ASDVBP60003008 ASDVBP6000508 ASDVBP6000308	AppliChrom ABOA StyDiViBe-P-10 <sup>6</sup> Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10.000Da - 4MioDa
ASDVBP70002508 ASDVBP70003008 ASDVBP7000508 ASDVBP7000308	AppliChrom ABOA StyDiViBe-P-10 <sup>7</sup> Å	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100.000Da - 10MioDa




Runge picture „homemade“

## AppliChrom ABOA CatPhil-P

Special GPC / SEC columns for the aqueous analysis of neutral, anionic and additionally cationic polymers.

For GPC / SEC analysis of

- Polycations, polyammonium compounds
- polyamines (chitosans), polyethylenimines,
- PEGylated polyethyleneimines
- polysaccharides
- polyanions (heparins, pectins, ...)




### Analyte: Poly DADMAC

Polydiallyldimethylammonium chloride, polyquaternium-6, Mw = 100-200kDa, CAS [26062-79-3]

Column: AppliChrom ABOA CatPhil-P-100  
AppliChrom ABOA CatPhil-P-350


Dimension: e.a. 300mm x 8mm  
Mobil Phase:  $\text{NaNO}_3$  + 0,2% formic acid in  $\text{H}_2\text{O}$   
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 100µl sample



### Analyte: High molecular weight chitosan, poliglusam, polyglucosamine, CAS [9012-76-4]

Column: AppliChrom ABOA CatPhil-P-400

Dimension: 3x 300mm x 8mm  
Mobil Phase:  $\text{NaNO}_3$  + 0,2% formic acid in  $\text{H}_2\text{O}$   
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 100µl sample




### Analyte: Chitosan sulfate

Column: AppliChrom ABOA CatPhil-P-100  
AppliChrom ABOA CatPhil-P-350

Dimension: e.a. 300mm x 8mm  
Mobil Phase: 0,05M  $\text{Na}_2\text{HPO}_4$  + 0,1M  $\text{NaNO}_3$  in  $\text{H}_2\text{O}$   
Flow: 1,0ml/min  
Temperature: 20°C  
Detection: RI  
Injection: 100µl sample

For aqueous polycations including

**Analyte: PEI (polyethyleneimine)**

linear, PEI 150, CAS [9002-98-6]

Mn = 4x10<sup>4</sup>, Mw = 4.2x10<sup>4</sup>,Mp = 3.9x10<sup>4</sup>Da;**PEI (polyethyleneimine)**

linear, PEI 25, CAS [9002-98-6]

Mn = 2.1x10<sup>3</sup>, Mw = 2.9x10<sup>3</sup>,Mp = 2.15x10<sup>3</sup> Da

Column: AppliChrom ABOA CatPhil-P-100  
AppliChrom ABOA CatPhil-P-350

Dimension: e.a. 300mm x 8mm

Mobil Phase: 0.1M NaCl + 0.2% TFA in H<sub>2</sub>O


Flow: 1.0ml/min

Temperature: 20°C

Detection: RI

Injection: 100µl sample

AppliChrom ABOA CatPhil-P Calibration curve  
Porosity P-100, P-200, P-350 and P400,  
based on dextrans

**Analyte: Dextrans**

Column: AppliChrom ABOA CatPhil-P-100  
AppliChrom ABOA CatPhil-P-200  
AppliChrom ABOA CatPhil-P-350  
AppliChrom ABOA CatPhil-P-400

Dimension: e.a. 300mm x 8mm


Mobil Phase: H<sub>2</sub>O,

Flow: 1.0ml/min

Temperature: 20°C

Detection: RI


Injection: 100µl sample



Polyethylenimine (PEI)


## AppliChrom ABOA CatPhil-P

Catalog #	Description	Dimension	Separation Range
SACP1002508 SACP1003008 SACP100508 SACP100308	AppliChrom ABOA CatPhil-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SACP2002508 SACP2003008 SACP200508 SACP200308	AppliChrom ABOA CatPhil-P-200	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-20.000Da
SACP3502508 SACP3503008 SACP350508 SACP350308	AppliChrom ABOA CatPhil-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	2.500Da-1.000.000Da
SACP4002508 SACP4003008 SACP400508 SACP400308	AppliChrom ABOA CatPhil-P-400	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10.000Da-5.000.000Da



## AppliChrom® ABOA CatPhil-P-JLJ series - additionally enhanced hydrophilicity

Special GPC / SEC columns for the aqueous analysis of neutral, anionic and additionally cationic polymers - JLJ series with further increased hydrophilicity for extended application range - also in pure aqueous eluents (calibration with dextran / pullulan and new: additionally with PEO / PEG or p-2-vinylpyridine possible).


**Analyte: Poly(2-vinylpyridine)**

CAS [25014-15-7], Mw = 40.000Da

Column: AppliChrom ABOA CatPhil-P-350-JLJ

Dimension: 300mm x 8mm


Mobil Phase: 0.1M NaCl + 0,2% TFA in H<sub>2</sub>O,

Flow: 1,0ml/min

Temperature: 20°C

Detection: RI

Injection: 20µl sample


**Analyte: Polyquaternium-33**

very high molecular weight, copolymer of trimethylaminoethyl acrylate salt and acrylamide, CAS [69418-26-4] (red) vs. high molecular weight polyacrylamide, CAS [9003-05-8] 15Mio Da (blue, dashed)

Column: AppliChrom ABOA CatPhil-P-100-JLJ

AppliChrom ABOA CatPhil-P-350-JLJ

AppliChrom ABOA CatPhil-P-500-JLJ

Dimension: e.a. 300mm x 8mm


Mobil Phase: 0.1M NaCl + 0,2% TFA in H<sub>2</sub>O,

Flow: 1,0ml/min

Temperature: 40°C

Detection: RI

Injection: 100µl sample


**Analyte: Poly DADMAC**

**Polydiallyldimethylammonium chloride, polyquaternium-6**

Mw = 100-200kDa, CAS [26062-79-3]

Column: AppliChrom ABOA CatPhil-P-100-JLJ

AppliChrom ABOA CatPhil-P-350-JLJ

AppliChrom ABOA CatPhil-P-500-JLJ

Dimension: e.a. 300mm x 8mm


Mobil Phase: 0.1M NaCl + 0,2% TFA in H<sub>2</sub>O,

Flow: 1,0ml/min

Temperature: 40°C

Detection: RI

Injection: 100µl sample




**Analyte:** Poly(vinylimidazole),  
CAS [25232-42-2]

**Column:** AppliChrom ABOA CatPhil-P-350-JLJ

**Dimension:** 300mm x 8mm  
**Mobil Phase:** 0.1M NaCl + 0,2% TFA in H<sub>2</sub>O,  
**Flow:** 1,0ml/min  
**Temperature:** 20°C  
**Detection:** RI  
**Injection:** 20µl sample

AppliChrom ABOA CatPhil-P-JLJ Calibration curve  
 Porosity P-100, P-200, P-350 and P400,  
 based on dextrans



**Analyte:** PEO  
 Dextran  
 Polyacrylamide

**Column:** AppliChrom ABOA CatPhil-P-100-JLJ  
 AppliChrom ABOA CatPhil-P-350-JLJ  
 AppliChrom ABOA CatPhil-P-500-JLJ

**Dimension:** e.a. 300mm x 8mm  
**Mobil Phase:** H<sub>2</sub>O,  
**Flow:** 1,0ml/min  
**Temperature:** 20°C  
**Detection:** RI  
**Injection:** 100µl sample

## AppliChrom ABOA CatPhil-P-JLJ

Catalog #	Description	Dimension	Separation Range
SACP100JLJ2508 SACP100JLJ3008 SACP100JLJ508 SACP100JLJ308	AppliChrom ABOA CatPhil-P-100-JLJ	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SACP350JLJ2508 SACP350JLJ3008 SACP350JLJ508 SACP350JLJ308	AppliChrom ABOA CatPhil-P-350-JLJ	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	2.500Da-1.000.000Da
SACP500JLJ2508 SACP500JLJ3008 SACP500JLJ508 SACP500JLJ308	AppliChrom ABOA CatPhil-P-500-JLJ	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10.000Da- ~50.000.000Da

A general guarantee for the functionality with all analytes can not be given. In individual cases it is always advisable to discuss the issue with our product specialists. An application laboratory in Oranienburg is also available for your service. Please contact us.

Errors and changes excepted. Products for research and development. Not tested for clinical, diagnostic or food applications.



## AppliChrom® ABOA HFIP-Phil-P

Special GPC / SEC columns for the GPC analysis of HFIP-soluble polymers.

For GPC / SEC analyzes from \*)

Polyesters (polybutylene terephthalate / PBT / CAS 24968-12-5,  
polyethylene terephthalate / PET / CAS 25038-59-9, polylactide  
PLA / CAS 26100-51-6)

Polyamide 6 / PA6 / polycaprolactam / CAS 25038-54-4, polyamide 6-6 /  
PA6-6 / polyhexamethylene adipamide / CAS131-17-2, polyamide 6-10 /  
PA6-10 / poly(hexamethylene sebacamide) / PA 6-10, CAS 9011-52-3

Other (paraformaldehyde / polyoxymethylene / POM / polyacetal / CAS  
30525-89-4 polyethylenimine / PEI / poly (iminoethylene / polyaziridine /  
CAS 9002-98-6)

HFIP GPC / SEC calibration versus

- PMMA / polymethylmethacrylate / CAS 901-14-7


HFIP Accessories:

5mM CF<sub>3</sub>COOK or 0.075M CF<sub>3</sub>COONa are added to the eluent HFIP  
for the suppression of electrostatic interactions and for artefact-free GPC.

Benefits AppliChrom ABOA HFIP-Phil-P GPC columns series over standard columns

- very large pore volume for high GPC resolution
- high GPC resolution for oligomers / condensates of 100-70,000Da
- high resolution GPC separation for the range 100-800.000Da
- large areas with high linearity calibratable areas without porosity artifacts
- low bleed
- standard columns also used for GPC-LS and GPC viscosity


\*) GPC / SEC of polylactide (PLA) - see also AppliChrom ABOA StyDiViBe in CHCl<sub>3</sub>, for polyethylenimine (PEI)  
- see also: AppliChrom ABOA CatPhil-P in H<sub>2</sub>O



**Analyte:** **Polymethylmethacrylat**  
 (PMMA), CAS 9011-14-7, CAS131-17-2  
 M = 901.000, 96.760, 32.500, 3196, 540Da

Column: AppliChrom ABOA HFIP-Phil-P-350


Dimension: 2x 300mm x 8mm  
 Mobil Phase: HFIP, 5mM CF<sub>3</sub>COONa  
 Flow: 0,5ml/min  
 Temperature: 40°C  
 Detection: RI  
 Injection: 100µl sample



**Analyte:** **Polyethylenterephthalat**  
 (PET), CAS 25038-59-9

Column: AppliChrom ABOA HFIP-Phil-P-350


Dimension: 2x 300mm x 8mm  
 Mobil Phase: HFIP, 5mM CF<sub>3</sub>COONa  
 Flow: 0,5ml/min  
 Temperature: 40°C  
 Detection: RI  
 Injection: 100µl sample



**Analyte:** **Polyethylenterephthalat**  
 (PET), CAS 25038-59-9

Column: AppliChrom ABOA HFIP-Phil-P-350

Dimension: 2x 300mm x 8mm  
 Mobil Phase: HFIP, 5mM CF<sub>3</sub>COONa  
 Flow: 0,5ml/min  
 Temperature: 40°C  
 Detection: RI  
 Injection: 100µl sample



**Analyte:** **Polyamide 6,6 (PA6-6),**  
 CAS131-17-2

Column: AppliChrom ABOA HFIP-Phil-P-350

Dimension: 2x 300mm x 8mm  
 Mobil Phase: HFIP, 5mM CF<sub>3</sub>COOK  
 Flow: 0,5ml/min  
 Temperature: 40°C  
 Detection: RI  
 Injection: 100µl sample; 1g/l

Good for many HFIP soluble polymers


## AppliChrom ABOA HFIP-Phil-P

Catalog #	Description	Dimension	Separation Range
SAHFIP1002508 SAHFIP1003008 SAHFIP100508 SAHFIP100308	AppliChrom ABOA HFIP-Phil-P-100	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	100Da-2.500Da
SAHFIP3502508 SAHFIP3503008 SAHFIP350508 SAHFIP350308	AppliChrom ABOA HFIP-Phil-P-350	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	up to 1.000.000 Da
SAHFIP5002508 SAHFIP5003008 SAHFIP500508 SAHFIP00308	AppliChrom ABOA HFIP-Phil-P-500	250mm x 8mm 300mm x 8mm 50mm x 8mm 30mm x 8mm	10.000Da- ~50.000.000Da

## Principle of GPC/SEC

### In detail:

Small polymers intrude to nearly all pores  
 Medium sized polymers intrude only to some pores  
 Very large polymers are totally excluded from the pores  
 => separation of polymers according to molecular size



**AppliChrom ABOA GPC/SEC Series:**

AppliChrom ABOA GPC columns are designed for high performance GPC separations. These columns contain high crosslinked polymeric particles of high reproducibility, good selectivity and high stability. Before shipment, each AppliChrom ABOA GPC column is tested and supplied with its individual test chromatogram and specifications. For long live time of AppliChrom ABOA GPC columns and for long term satisfaction it is important for customers to take care about some points that will be given in the following text. Due to several different specific applications of customers, the text can not be absolute – therefore, if you have any question, please do not hesitate asking the AppliChrom service at [sales@applichrom.de](mailto:sales@applichrom.de) and inform your local distributor.

After receiving of column:

- Control if there is any obvious physical damage from shipping
- Check if the column you got is the column you ordered
- All columns are shipped in the solvent that is documented in the column test certificate
- Verify the column quality

**Please take into account that performance of column also depends on the chromatographic system at customers site.**

**Following parameters can have an influence the verification of column performance:**

- Loop size, degree of loop filling, concentration of sample
- Capillary diameters and length in chromatographic system
- Temperature of column
- Effective detector cell volume
- Data acquisition rate
- Fitting connections

**AppliChrom ABOA GPC column series connectors:**

- In case a set of columns of different porosities plus eventual precolumn is ordered, customer will be also supplied with the suited column connectors, of 1/16" OD. The fittings have the 10-32 type thread and ferrules that are pre-fixed to ensure tightening without voidvolume

**Tubing and fittings:**

- 8mm ID AppliChrom ABOA GPC column series are designed to be connected with 1/16" OD capillary via 10-32er threading with HPLC-System
- If connecting chromatographic system with columns, please take care about the depth of capillary coming from system - connected by ferrule and fitting to columns. If distance between capillary end and ferrule is too low you get extra void volumes and loss of peak performance, if distance between capillary and ferrule is too long, system can leak or column inlet can be destroyed
- **See also: AppliChrom HPLC / GPC column user guide**

## Installation of column

- Rinse chromatographic system (without column) with degassed and filtered eluent (here: DMSO eluent) with a moderate flowrate, check if there is no leak in the system, there is no blockade in the system and please make sure that detector will not be blocked or getting a too high pressure! Maybe laboratory must be continuous good temperatured (warm enough) to be sure DMSO is still fluid and viscosity is low enough! Please ensure this temperature also at weekends and/or holidays. Please document system pressure for the later used flowrate, at e.g. 0.4ml/min. If salts are added to the eluent, please take care to avoid precipitation of salt in system – as a result of this irreversible damages can occur
- Connect column respective columns with capillary from injector. Take care columns are connected in flow direction. The outlet of column (of the final column if there is a set of column) should be connected direct with waste, bypassing the detector system
- If using more than 1 column (series of columns); take the columns in the following order:  
inlet => precolumn, maincolumn of smallest porosity, maincolumn of medium porosity, maincolumn of largest porosity => outlet
- Heat the columns at least to 50-60°C or to about 80°C for at least 1h
- Start pumping with 0.1ml/min, check if connections are tight, observe the outlet of columns to make sure a clear eluent flows out. Pump at least 5min / column at this flowrate, document the observed pressure.
- Increase stepwise to the final flowrate, e.g. 0.4ml/min, document the backpressure
- If backpressure is constant and eluent leaving the outlet pipe is clear, connect with detector (please make sure that connector is warm enough for ensuring fluidity of DMSO)
- After a total approx 5-10 column volumes baseline should be stable and measurements can start.
- If system is not needed for some time (over night, weekends,...), please lower flowrate e.g. to 0.05 or 0.1ml min, but do not change the column temperature. Every cooling and heating process of column means stress for column
- Please make sure that there is not a column pressure drop of more than 80bar per column or 30bar per Precolumn
- Avoid fast pressure increase or fast pressure lowering of columns
- If you inject sample, please take care sample is clear and free of particles (e.g. by filtration) and please take the lowest concentration that is necessary to get the results you need. Especially if using / analysing reactive samples, column lifetime can be increased significantly by using low concentrations of sample. Samples should be diluted with eluent to avoid eventual extra-peaks

## Important:

**If you have any question, please do not hesitate contacting the AppliChrom service, e.g. at [sales@applichrom.de](mailto:sales@applichrom.de)**

**We want to get you as satisfied customer with good products plus good service.**

## Guarantee:

AppliChrom guarantees the quality according to its test certificate at the date of shipment. Column warranty is voided if end-fittings are removed. Applications at customer's site are not in the responsibility of supplier.

All information on these sides are put together with the biggest care, has been worked on and been checked. No liability or guarantee can be still assumed for content and objective correctness as well as actuality and completeness. All information without guarantee.


**Installation:**

- a) Rinse first your HPLC / GPC system (suction tube, pump, injector incl. injection loop as well as capillary connections up to the column entrance) with the mobile phase, e.g., with 1ml / min, note the originating counterpressure.
- b) Connect the capillary to the column entrance with the HPLC / GPC column. Provide for the fact that the capillary is optimally connected with the column entrance. Between the capillary end and the column should be
  - no dead volume originate (by too short capillary piece) – leads to line widening
  - the capillary end also not longer as the admission within the HPLC / GPC-Column – can lead to leakage in the column head.

**Please follow,**

the space of the capillary ends can vary between different HPLC / GPC suppliers. To make sure that the correct space of the capillary end are reached, the adaptation of the ferrule is recommended in respective HPLC / GPC hardware.

Charting of the influence of the capillary length



- c) **Pump** the eluent by the HPLC / GPC system consisting of pump, injector, capillary between column and injector, column, e.g., with 1ml / min with 4.6 mm of columns ID. Let go the first eluent 2-4ml then the column leave directly in the waste - under avoidance of the detectors (is valid for columns to max. 4.6 mm ID and you let go lengths from 50-250mm), with 8 mm of columns ID and lengths from 250-300mm please 15-20ml eluent directly in the waste – under avoidance of the detectors. If you control and document, besides, please pressure, rinse the column in this manner so long to the measured pressure is steady and the eluent clearly of the column flows out.
- d) **Connect** the column end from that of the eluent resigned with the detector system. Guarantee that yourselves, besides, no significantly raised system counterpressure compared with the before measured value from c) puts (tip to blockage of the detector system, e.g., from dried up eluent, tough or crystallised solvents, mature test, blocked, narrowed or broken off capillaries...) - blockage of the detector system can lead to the destruction with pressure excess – use please only accordingly trained, attentive staff!
- e) **Use** the fitting and ferrule connection between column and detector system in the manner that a) **no dead volume** (line widening) and b) no excess end of the inside-recumbent capillary originates (can lead to the destruction of the column exit). Ideally is the new touch-down of a connection from ferrule and fitting or, e.g., the use of PEEK-screw connections at the column exit.
- f) **Document** at regular intervals the system pressure.

**Please, use for AppliChrom-precolumnd cartridges only the AppliChrom-precolumnd cartridge holder fitting for it!**

All information on these sides are put together with the biggest care, has been worked on and been checked. No liability or guarantee can be still assumed for content and objective correctness as well as actuality and completeness. All information without guarantee.

## General terms and conditions

The prices are free, in EURO, without VAT, this is calculated separately.

The listed prices are valid in the Federal Republic of Germany, deviations may result from changes in the manufacturers' prices.

Our offer is a business-to-business (B2B) offer and also addressed to selected traders and freelancers as well as universities and other public clients but not to private customers.

### Delivery terms:

Federal Republic of Germany: Delivery is not free (EXW Oranienburg, INCOTERMS 2010). We charge a shipping and packaging fee depending on the size and weight of the goods or by arrangement.

(EU and third countries): Delivery is carried out DDU / EXW Oranienburg, INCOTERMS 2010. If no specification of the buyer is made, we determine the freight carrier and charge the buyer the shipping costs as well as all additionally incurred costs.

### Payment:

From invoice date within 30 days net or by appointment.

### Complaints:

Are possible within 8 business days after receipt of the goods.

In any case, they must be notified in writing.

The return of the goods is only permitted with our consent.

In case of justified defects, we have the choice between exchange of the goods, rework or reimbursement of the purchase price.

Complaints do not entitle the buyer to suspend his payment.

Any further damages claims of the buyer are excluded, as long as they are not caused by gross negligence or intentional. In case of gross negligence or intentional breach of our contractual or statutory obligations, our liability is limited to the foreseeable damage; we do not assume any liability for any defects.

Claims under the Product Liability Act are not affected by the above regulation.

If the price offer is obviously based on an incorrect price due to transmission errors or error, we are entitled to make a subsequent price correction.

If the customer rejects this correction or if no agreement can be reached between the parties to the contract, we are entitled to a right of withdrawal at any time, as far as legally permissible.

### Warranty:

Corresponds to legal requirements.

However, the guarantee / warranty of the manufacturer applies to commercial products.

The guarantee for the properties of the columns is limited to the conditions of the original test chromatogram.

AppliChrom ensures the quality of the AppliChrom HPLC / GPC columns under the following conditions:

- If a column does not function properly, AppliChrom will re-install the column or replace it without incurring any costs for the customer.
- In the case of return of columns, AppliChrom must first issue a return authorization. Defective products must be accompanied by a written declaration of the defect. A redemption only takes place under the following conditions:
  - All columns must be checked upon receipt and AppliChrom must have all defects within 14 days after receipt of the column.
- Unless otherwise agreed, the warranty period is limited to a maximum of 90 days after delivery date.
- The guarantee for the properties of the columns is limited to the conditions of the original test chromatogram.

### To void the warranty / guarantee:

- Removing or unscrewing the end fittings of the HPLC/GPC columns automatically voids column warranty.
- Physical damage to the column due to misuse or mishandling, including mechanical shocks.
- Damage to the stationary phase or the base material due to handling in the case of incorrect chemical conditions, unsuited solvent, temperatures or pressure conditions.
- Failure due to high backpressures caused by improper solvent and/or sample filtration practices causing particulate build-up or precipitation in the column or end-fitting.
- Incorrect selection of packing material made by customer for their particular use or incompatibility of equipment, etc.

### Use restriction

Only for R & D as well as laboratory, not tested for pharmacy or medical diagnostics.

Buyer will not use AppliChrom's name in any publicity or advertising without AppliChrom's prior written approval.

Proprietary right of reservation exists until full payment of the goods.

Only place of performance and jurisdiction is Oranienburg/Germany.

We reserve the right to make technical changes.

Our company does not participate in a consignment procedure.



Application	Page	Application	Page
Acrylate water-based (range: 100-1.000.000Da)	14	PET, Polyethylene terephthalate (in HFIP)	40
Alginate-Na (range: 100-1.000.000Da)	15	PMMA, Polymethylmethacrylate (in HFIP)	40
Aliphatic polyester	30	PMMA, Polymethylmethacrylate (in THF)	31
Bisphenol-A-Epichlorohydrin (I) (pore perfect)	29	PMMA, Polymethylmethacrylate, oligomer	28
Bisphenol-A-Epichlorohydrin (I) (pore to large)	29	Poly DADMAC	34
Bisphenol-A-Epichlorohydrin (I) (pore to small)	29	Poly DADMAC (comparison 3 samples)	37
Caramel color GPC (range: 100-1.500.000Da)	10	Poly(2-vinylpyridine), Mw = 40.000Da	37
Carrageenan (range: 100-1.000.000Da)	15	Poly(N-isopropylacrylamide)	7
Casein fraction water-soluble	17	Poly(vinylimidazole)	38
Chitosan sulfate	34	Poly[2-(butenyl)2-oxazoline-co-ethylenimine]	8
Chitosan, high molecular weight	34	Polyacrylamide, comparison of PEO, dextrans	38
Comparison of 3 samples, Poly DADMAC	37	Polyamide 6,6, (PA6-6) (in HFIP)	40
Comparison of PEO, Dextrans, Polyacrylamide	38	Polybutyl methacrylate/	
Corn syrup (range: 100-1.000.000Da)	14	Poly(dimethylamino-ethylmethacrylate) M=22kDa	9
Dextran 650	6	Polycarboxylate ether	16
Dextran sulfate-Na (range: 100-1.000.000Da)	15	Polydiallyldimethylammonium chloride,	
Dextrans	35	polyquaternium-6,	34
Dextrans, comparison of PEO, Polyacrylamide	38	Polyester, aliphatic	30
DP2, mixture of...4	21	Polyether polyol	30
DP3, mixture of...4	21	Polyethylene glycol (PEG)	28
Ethanol, mixture of...3	21	Polyethylene terephthalate, PET (in HFIP)	40
Fish protein hydrolysates (range: 100-70.000Da)	16	Polyethyleneimine, PEI	35
Fructose, mixture of...1*	21	Polyglucosamine, high molecular weight	34
Fructose, mixture of...2*	21	Polyglusam, high molecular weight	34
Fructose, mixture of...3*	21	Polymethylmethacrylate, PMMA (in HFIP)	40
Fructose, mixture of...4*	21	Polyol soft foam	30
Glucose, mixture of...1*	21	Polyquaternium-33	37
Glucose, mixture of...2*	21	Polysaccharide	6
Glucose, mixture of...3*	21	Polystyrene (PS), Mp = 578Da	28
Glucose, mixture of...4*	21	Polystyrene, PS (in THF)	32
Glycerin, mixture of...1*	21	Polyvinyl alcohol M=22kDa	8
Glycerin, mixture of...2*	21	Polyvinyl alcohols, 88% degree of hydrolysis	17
Glycerin, mixture of...3*	21	Polyvinyl pyridine degraded	8
Grapeoil (Triglyceride GPC)	28	Polyvinyl pyridine fractions	8
Hard Foam, polyol	30	Polyvinylchloride, PVC (in THF)	31
Heparin.Na, 8-25kDa	15	Pork gelatin (range: 100-1.000.000Da)	18
Honey	22	Pork gelatin vs. Gelatin from collagen hydrolysate	18
Humic acids / humates (range: 100-1.500.000Da)	10	PS, Polystyrene (in THF)	32
Hyaluronic acid	14	PS, Polystyrene, Mp = 578Da	28
Inulin (range: 100-1.000.000Da)	14	Pullulan	16
Jelly Bean „Gumminbärchen“	17	PVC, polyvinylchloride (in THF)	31
Lignin conifer bark	6	PVP, Polyvinyl pyrrolidone	17
Malt Beer	22	Silicone (in Toluene)	31
Maltodextrin 6	25	Soft foam, polyol	30
Maltodextrin 6, detailed view	25	Sorbitol, mixture of...1	21
Maltodextrin 12	24	Soy protein GPC	
Maltodextrin 12, detailed view	24	(range:100-1.00.00Da (based on dextrans))	9
Maltodextrin 12 (range: 100-1.000.000Da)	16	Spruce bark extract	6
Maltodextrin 19	24	Starch hydrolysates	13
Maltodextrin 19, detailed view	24	Sucrose (saccharose), mixture of...1	21
Mannitol, mixture of...1	21	Sucrose (saccharose), mixture of...2	21
Manuka honey protein GPC		Sucrose (saccharose), mixture of...3	21
(range 100-1.000.000Da (based on dextrans))	9	Triglyceride GPC, Grapeoil	28
Melamin (range: 100-1.500.000Da)	10	UF-resin	7
MUF-resin	7	Water-based acrylate (range: 100-1.000.000Da)	14
Oligosaccharide	13	Water-based, acrylate (range: 100-1.000.000Da)	14
PA6-6, Polyamide 6,6 (in HFIP)	40	Water-soluble casein fraction	17
Pea protein GPC		Water-soluble, Casein fraction	17
(range: 100-1.000.000Da (based on dextrans))	9		
Pea starch	7		
Pectin	13		
PEG / PEO	12		
PEG, Polyethylene glycol	28		
PEGylated protein	13		
PEI, polyethyleneimine	35		
PEO / PEG	12		
PEO, Dextrans, Polyacrylamide, comparison of	38		

\*

Mixture 1	Mixture 2	Mixture 3	Mixture 4
1. Sucrose	1. Sucrose	1. Sucrose	1. Dp 3
2. Glucose	2. Glucose	2. Glucose	2. Dp 2
3. Fructose	3. Fructose	3. Fructose	3. Glucose
4. Glycerin	4. Glycerin	4. Glycerin	4. Fructose
5. Mannitol		5. Ethanol	
6. Sorbitol			



